AT02004, AT02005, AT02006, AT02007

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Thank you for purchasing our product. Please check the condition of the device before doing anything. If the product operates or looks abnormal, immediately contact the dealer from whom you purchased the device. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions or instructions, contact a health care professional or your dealer to avoid damaging the product.

NOTE: Inspect all parts for damage caused in transit. If you notice such damage, DO NOT use the product. Contact the manufacturer for more information.

DESCRIPTION AND PURPOSE

Rehabilitation deambulatorettes are lightweight, comfortable-to-use products designed to help people with limited mobility get around. They provide the user with support, reducing the risk of falling over, and allow walking with support on both sides.

MODELS:

Depending on the needs and preferences of the user, it is possible to choose from 4 versions of rehabilitation rollators:

AT02004 - version with 4 feet, adjustable. It provides the greatest stability. In order to rearrange the rehabilitation rollator, it is necessary to raise it. It is mainly recommended for moving short distances.

AT02005 - version with 2 wheels. Allows you to move the rehabilitation rollator without lifting it.

AT02006 - version with 4 wheels. Rear swivel wheels with brakes. The product allows the user to move forward quickly and, if necessary, it is possible to lock the rehabilitation rollator by using the brakes. The model is recommended for people with no balance problems.

AT02007 - version with 4 swivel wheels. The product allows the user to move quickly in all directions, and if necessary, it is possible to block the rehabilitation rollator by activating the brakes. This model is recommended for people with high mobility and no balance problems.

CONTRAINDICATIONS

physical or mental limitations (e.g., visual impairment) that prevent safe handling.

1. severe brain injuries: In cases of severe brain injury, especially when there is significant loss of motor function, a rehabilitation walker may be inappropriate or even dangerous.
2. Significant stability problems: People with significant body stability problems may have difficulty controlling a rehabilitation walker and may be more likely to fall or be injured.
3. Balance incontinence: If a person has difficulty maintaining balance or is prone to falls, using a rehabilitation walker can be dangerous.
4. Coordination problems: People with severe coordination problems may find it difficult to use a rehabilitation walker effectively and may not get the expected benefits from its use.
5. Postural defects: People with postural defects, such as scoliosis, may require special assessment and adaptation of the rehabilitation walker to avoid worsening their condition.

WARNINGS:

- 1) Do not use the AT02006 ora with the AT02007 for people with balance disorders. This may lead to the user falling over.
- 2) Do not lean outside the perimeter of the rehabilitation rollator. This risks the user falling over.
- 3) Do not rock the rehabilitation rollator. This risks the user falling over.
- 4) In the case of a four-wheeled rehabilitation rollator (AT02006, AT02007), when not in use at any given time - its brakes must be engaged.
- 5) Do not use if the weight of the user exceeds 150 kg.
- 6) Avoid contact with open flame. The product may catch fire.
- 7) Make sure that after installation, the steel pin of the foot or caster mounting clip passes through the hole on both sides of the wall. Otherwise, the product may not be stable and the hole may become deformed. Also, the pin that changes the rehabilitation rollator from walking to fixed should pass through the hole on both sides of the wall.

ATTENTION: it is forbidden to use the product other than for its intended purpose.

NOTE: Before use, check that all components have been assembled correctly.

NOTE: The manufacturer is not responsible for damage caused by negligent maintenance, inadequate servicing or resulting from failure to follow the recommendations in this manual.

CAUTION: If the product is not used properly, there may be a risk of tipping over. Please follow the instructions for moving.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this device is intended (see the section on the purpose of the device of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device independently or as indicated by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer, should be taken into account.

NOTES:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person or
- (b) temporary or permanent deterioration of the health of the patient, user or other person or
- (c) a serious threat to public health


the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions, or other disturbing symptoms that are unclear to you related to the use of the medical device, consult a healthcare professional.

PREPARATION FOR USE:

The rehabilitation rollator is delivered folded for transport.

DISCLAIMER

The rehabilitation rollator is unfolded by grabbing the front part of the frame near both cord attachment points and pulling the two parts away from each other until they stop.

SUBMITTING

To fold the rehabilitation rollator, pull the ball (knob) inward.

HEIGHT ADJUSTMENT

Before using the product, adjust its height. To do this, remove the clip, adjust the height - the same for all four legs, and then slide the clip into the corresponding hole. **MAKE SURE THAT THE STEEL PIN GOES THROUGH THE HOLE ON BOTH SIDES.** Otherwise, the product may not be stable, and the hole may be deformed.

MOVING AROUND WITH A REHABILITATION ROLLATOR

AT02004 - The rehabilitation rollator is repositioned by lifting it and placing it slightly in front.

AT02005 - The rehabilitation rollator moves forward by lifting its back, pushing it and placing the back part on the ground.

AT02006 - The rehabilitation rollator is repositioned in any direction by pushing it. IF THE USER IS NOT USING THE PRODUCT AT THE MOMENT, THE BRAKES LOCATED AT THE WHEELS SHOULD BE ACTIVATED (BY PRESSING).

AT02007 - The rehabilitation rollator is repositioned in any direction by pushing it. IF THE USER IS NOT USING THE PRODUCT AT THE MOMENT, THE BRAKES LOCATED AT THE WHEELS SHOULD BE ACTIVATED (BY PRESSING).

This symbol indicates the maximum weight of the user.

SPECIFICATION AT02004

Width (front): 47 cm

Width at widest point: 58 cm

Width when folded: 25 cm

Width between handles: 41 cm

Height: 80-93 cm

Length: 60 cm

Weight: 1.4 kg

Maximum user weight: 150 kg

SPECIFICATION AT02005

Width (front): 47 cm

Width at widest point: 58 cm

Width when folded: 25 cm

Width between handles: 41 cm

Height: 80-93 cm

Length: 60 cm

Weight: 1.6 kg

Maximum user weight: 150 kg

SPECIFICATION AT02006

Width (front): 50 cm
Width at widest point: 58 cm
Width when folded: 25 cm
Width between handles: 41 cm
Height: 80-93 cm
Length: 60 cm
Weight: 2.3 kg
Maximum user weight: 150 kg

SPECIFICATION AT02007

Width (front): 50 cm
Width at widest point: 58 cm
Width when folded: 25 cm
Width between handles: 41 cm
Height: 80-93 cm
Length: 60 cm
Weight: 2.75 kg
Maximum user weight: 150 kg

COMPONENTS: Rehabilitation rollator, instruction manual.

MAINTENANCE AND CLEANING OF THE PRODUCT

The product should be cleaned periodically with water and cleaning agents. The product can be treated with disinfectants. Do not use rough objects for cleaning.

STORAGE

When not in use, the product should be stored in a dry and cool room away from direct sunlight. Avoid prolonged exposure to direct sunlight.

DISPOSAL OF THE PRODUCT

Once the product is taken out of service, it can be treated as municipal waste.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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