AT03551/AT03552(S)/AT03553
INSTRUKCJA UŻYTKOWANIA, KLIN REHABILITACYJNY
USER MANUAL , REHABILITATION WEDGE
GEBRAUCHSANWEISUNG, LAGERUNGSKEIL
NÁVOD K POUŽITÍ, REHABILITAČNÍ KLÍN
NÁVOD NA POUŽÍVANIE, REHABILITAČNÝ KLIN
MODE D'EMPLOI, CALE DE RÉHABILITATION
GEBRUIKSAANWIJZING, REVALIDATIEKUSSEN
INSTRUCCIONES DE USO, CUÑA DE REHABILITACIÓN
ISTRUZIONI PER L'USO, CUEO RIABILITATIVO
BRUKSANVISNING, REHABILITERINGSKIL

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EN

Thank you for purchasing our product. Please read the user manual carefully before use. Do not use the product until you have read and understood this manual. If you do not understand the warnings, notes, and recommendations, contact a healthcare specialist or the retailer to avoid product damage or personal injury.

WARNING: Check all parts for damage incurred during transportation. If such damage is detected, DO NOT use the product. Contact the retailer for further information.

TARGET PATIENT GROUP

Individuals dealing with diseases, dysfunctions, or injuries for treatment/rehabilitation or compensation for which this product is intended (see the section on the product's purpose in this user manual). The product can be purchased independently by the user or based on recommendations from a doctor, therapist, or other specialist. Whether purchasing the product independently or following a recommendation from a doctor/therapist/other specialist, it is essential to consider the available sizes/necessary functions and product variants, indications, and contraindications for use, as well as the information provided by the manufacturer.

APPLICATIONS/INDICATIONS

The rehabilitation wedge is a tool used in therapy and health prevention to support proper body positioning, stabilization, and relief for specific muscle and joint areas. Its use is particularly popular in rehabilitation, physiotherapy, and daily care for immobilized individuals.

Applications of the rehabilitation wedge:

Stabilization and support of the body:

The rehabilitation wedge helps stabilize and appropriately position different parts of the body during exercises or rest.

It provides support for the lumbar spine, knees, ankles, hips, and other areas requiring relief.

Improvement of body posture:

Used to correct body posture, particularly in situations where the patient requires proper alignment, e.g., in lying or sitting positions.

Supports the correct positioning of the spine and pelvis, which helps prevent back pain and postural problems.

Rehabilitation after surgeries and injuries:

The wedge is useful in the rehabilitation process following orthopedic surgeries, such as hip or knee joint replacements.

It helps relieve stress on damaged areas and maintain proper body positioning during rehabilitation exercises.

Reduction of muscle tension:

The rehabilitation wedge helps relax tense muscles, especially in the lower back, neck, and shoulders, by providing appropriate support and body positioning.

It can be used during stretching or relaxation exercises to increase comfort and effectiveness.

Facilitation of positions for rehabilitation exercises:

The wedge is helpful during various rehabilitation exercises, especially those requiring elevation or stabilization of specific body parts.

It assists in controlled stretching and muscle strengthening, supporting the process of regaining full functionality.

Indications for the use of the rehabilitation wedge:

Pain in the spine and joints:

Recommended for individuals suffering from chronic back pain, especially in the lumbar and cervical regions, related to improper posture.

Used for hip and knee joint pain to provide relief and support for these areas.

After orthopedic and neurological surgeries:

Used after surgical procedures, such as joint surgeries, to ensure proper body positioning and reduce tension in operated areas.

Helps patients after neurological injuries who require special care for proper body and muscle positioning.

Bedsores prevention:

The rehabilitation wedge supports bedsores prevention in immobilized individuals or those with limited mobility by reducing pressure on sensitive body areas.

It is particularly recommended for bedridden patients, e.g., in a hospital or at home.

Problems with blood and lymph circulation:

Recommended for individuals with circulation issues, such as varicose veins, swelling, or poor circulation in the lower limbs.

Using the wedge under the legs supports blood and lymph flow, reducing the risk of swelling and discomfort.

Support during pregnancy:

The rehabilitation wedge is helpful for pregnant women, providing support for the back and legs, which helps alleviate pain and tension caused by the growing belly.

It can also be used postpartum to support muscle and joint recovery.

CONTRAINDICATIONS

Acute inflammation of muscles, joints, or soft tissues.

Recent injuries or open wounds in areas of contact with the wedge.

Blood circulation problems, such as thrombosis, which may require consultation with a doctor before use.

Caution, if the product is used by a person who is unable to move independently, it must be ensured that the user's body position relative to the product is regularly changed – otherwise, pressure sores may develop.

SPECIFICATION

-polyurethane foam

-foam does not deform

-leather-like fabric cover

-easy to clean

AT03551 – 44 x 33 x 22 cm

AT03552 – 48 x 38 x 25 cm

AT03552(S) – 40 x 43 x 27 cm

AT03553 – 53 x 37 x 29 cm

USAGE INSTRUCTIONS

Preparation for use:

Choose a flat, stable surface for using the wedge, such as a bed, exercise mat, or couch.

Ensure the wedge is positioned correctly to provide comfort and stability during use.

Body positioning on the wedge:

Support for the lower back or hips: Lie on your back and place the wedge under your lower spine or hips. This will gently elevate the hips, reducing pressure on the lumbar spine and improving lying comfort.

Leg elevation: To support blood circulation in the legs, place the wedge under the knees or ankles, lifting the legs upward. This position helps reduce swelling and supports lymphatic drainage.

Neck and cervical spine support: You can also use the wedge under your neck and upper back to relieve cervical spine pressure and reduce neck muscle tension.

Positions for rehabilitation exercises:

When performing stretching or rehabilitation exercises, position the wedge to provide stable support for the body part in question, e.g., under the feet during leg strengthening exercises.

The wedge can also be used in a side-lying or seated position to provide the appropriate body angle during exercises.

Comfort and pressure control:

Ensure you feel comfortable in every position on the wedge. If you experience any discomfort, adjust the wedge's placement or change the angle.

There should be no feeling of pain or excessive pressure on any part of the body. If necessary, reduce the time spent in the position.

Duration of use:

Initially, use the wedge for short periods, about 10-15 minutes, to allow the body to adapt to the new position.

As your body adjusts to the wedge, you can gradually extend the duration of its use, depending on your therapist's recommendations or your own feelings.

Ending usage:

After completing a session with the wedge, slowly change position and get up from the mat or bed to avoid sudden muscle and joint strain.

Perform a few gentle stretches to restore the body's natural position and increase blood flow.

Important tips:

Regularity: For best results, use the wedge regularly according to your physiotherapist's recommendations.

Body reactions: Always monitor how your body responds to the wedge positions. If you experience pain or discomfort, stop using it and consult a specialist.

Hygiene: Keep the wedge clean by regularly washing it according to the manufacturer's instructions.

LIST OF COMPONENTS

rehabilitation wedge, user manual

WARNING:

In the event of a "serious incident" related to the product that directly or indirectly resulted in, could have resulted in, or may result in any of the following events:

a) death of the patient, user, or another person, or

b) temporary or permanent deterioration of health of the patient, user, or another person, or

c) serious threat to public health,

the "serious incident" must be reported to the manufacturer and the appropriate authority of the member state where the user or patient resides. In Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

WARNING: If you experience pain, allergic reactions, or other worrisome or unclear symptoms related to the use of the medical device, consult a healthcare professional.

WARNING: Using the product in a manner inconsistent with its intended purpose is prohibited.

CLEANING

Cleaning the foam requires a gentle approach to avoid damaging its structure. Here is a guide on how to clean it properly:

What you will need:

Mild detergent or dishwashing liquid

Warm water

Baking soda (optional for deodorization)

Vacuum cleaner (preferably with a soft brush)

Spray bottle

Soft cloth or sponge

Enzymatic cleaner (optional for stain removal)

Cleaning steps:

Vacuuming the foam:

Remove dust and dirt: Use a vacuum cleaner with a soft brush to gently clean the foam surface. This removes dust, dirt, and debris before wet cleaning.

Spot cleaning with mild detergent:

Mix detergent: Combine a few drops of mild detergent with warm water in a spray bottle.

Lightly spray the foam: Avoid soaking the foam. Lightly spray the soiled area.

Blot the stains: Use a soft cloth or sponge to gently blot the area. Do not rub, as this may damage the foam's structure.

Rinse the area: If necessary, lightly spray with clean water to remove detergent residue.

Handling tough stains:

Biological stains or urine: Use an enzymatic cleaner. Apply it lightly to the stain, leave it for a moment (as per product instructions), then blot with a clean cloth.

Grease stains: Use a small amount of dishwashing liquid or a special upholstery cleaner designed for foam materials.

Deodorization (optional):

Use baking soda: Sprinkle a layer of baking soda on the foam and leave it for several hours or overnight to absorb odors.

Vacuum the baking soda after sufficient time.

Drying the foam:

Air dry completely: Foam is highly absorbent, so it may take some time to dry. Leave it in a well-ventilated area.

Use a fan to speed up the process: Never use a direct heat source, such as a hairdryer, as it may damage the foam. A cool fan is best.

Ensure the foam is completely dry before reuse to prevent mold or mildew growth.

What to avoid:

Never submerge the foam in water, as it may be difficult to dry and could become damaged.

Avoid using strong chemicals such as bleach or alcohol-based cleaners, which may weaken the foam's structure.

Avoid direct heat: Heat can damage the foam structure, so do not tumble dry or use hot air.

Cleaning the cover: Wash at 95°C. Do not bleach. Do not dry clean. Do not iron. Do not tumble dry.

STORAGE

Store the product in a dry place, away from direct sunlight. Not suitable for children under three years of age. Keep away from fire.

DISPOSAL

After decommissioning, the medical device can be disposed of as regular municipal waste. For proper disposal, please contact your local authority or the nearest waste collection point to learn how to dispose of the product correctly. By adhering to local regulations, you help protect the environment.

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Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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