AT04001, AT04002, AT04004

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

Thank you for purchasing our product, please read the instruction manual carefully. Do not use the product before reading and understanding these instructions. If you do not understand the warnings, cautions or instructions, please contact your healthcare professional or retailer to avoid damage to the product.

The Dessault shoulder and arm brace is made of lightweight, breathable foam trimmed with skin-friendly fabric. The material wicks heat and moisture to the outside, which increases the comfort of the product and makes it easy to maintain hygiene.

INTENDED USE

The brace is recommended for the treatment of injuries to the shoulder, elbow and shoulder joint, including fractures in the shoulder girdle, dislocations and sprains of the shoulder and elbow joint. It can also be used to treat inflammation and degenerative changes of the shoulder and elbow joint. The brace minimises the period of plaster immobilisation, guaranteeing perfect stabilisation and immobilisation of the shoulder joint and upper limb. The product has a wide range of adjustment. The brace can be quickly and easily adjusted to the patient's body shape. The Velcro fastening guarantees effective immobilisation and pressure relief of the upper limb. The brace should be used as recommended by a doctor or physiotherapist. The brace is reversible and can be used on the right or left arm. The product should not be used in case of open wounds, abrasions of the epidermis or after the previous application of warming gel or ointment.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/sizes/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

CAUTION:

In the event of a device-related "serious incident" that directly or indirectly led to, could have led to or may lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION:

If you experience pain, allergic reactions or other distressing symptoms that are unclear to you related to the use of the medical device, consult a healthcare professional.

ATTENTION: It is forbidden to use the product other than for its intended purpose.

CLEANING

The orthosis may only be washed by hand in lukewarm water (up to 40°C) with a mild detergent (e.g. grey soap). Do not iron, do not dry clean, do not use bleach, do not spin dry, do not tumble dry. Dry unfolded, squeeze gently before drying.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

DELIVERY: Orthosis, instructions for use

DISPOSAL OF USED PRODUCT


After withdrawal from use, the product may be treated as municipal waste.

SIZES

To choose the right size, measure your chest circumference:

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




WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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