AT51006

INSTRUKCJA OBSŁUGI CHODZIK
USER MANUAL ROLLATOR
BENUTZERHANDBUCH ROLLATOR
NÁVOD K POUŽITÍ CHODÍTKO
NÁVOD NA POUŽITIE CHODÍTKO
MANUEL DE L'UTILISATEUR ROLLATEUR
HANDLEIDING ROLLATOR
MANUAL DE USUARIO ANDADOR
MANUALE D'USO ROLLATORE
ANVÄNDARMANUAL ROLLATOR

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Thank you for purchasing our product. Before using, please read this user manual carefully. Incorrect use may lead to health detriment, injury or accidents. Therefore, in order to use the product safely and successfully, please read these instructions carefully.

NOTE: Inspect all parts for damage caused in transport. If any such damage is found, DO NOT use the product. More information from the manufacturer.

GENERAL WARNING

Do NOT assemble and use the product until you have read and understood this manual. If you do not understand any warnings, cautions, or instructions, contact your healthcare professional, distributor, or service technician before attempting to assemble the product to avoid personal injury or damage to the product. Each user should ask his doctor or therapist to determine the appropriate variant of regulation and method of use.

The rollator is not intended to be guided by the user while sitting on it. A doctor / therapist should be present when adjusting the height of a rollator to provide support for the user and to correctly adjust the brakes and tighten the screws. When the rollator is unfolded, always make sure that the seat tube is locked and properly attached to the frame, and that the brake cable is inside the side tube and is securely attached to it. Check the brake cable after each height adjustment to avoid loosening the cable. Make sure the knobs are locked so that the main support structure will not move as a result of pressure. The quick release mechanism is installed in place and will remain in the release position even if you pull forcefully in the opposite direction. The copper pins on the adjustment tube pop out from the right.

Failure to follow the above recommendations may lead to injuries. Before use, make sure that ALL the adjustment elements of the handles and the seat are secure and that the wheels and moving parts of the product are in good condition. During use, ALL the wheels must be in contact with the ground AT ALL TIMES to ensure the stability of the rollator. BEFORE the user is seated on the seat, the brakes MUST be in the locked position. When the rollator is used in a stationary position, its handbrakes MUST be in the locked position.

ATTENTION: The manufacturer is not responsible for damage caused by neglect of maintenance, inadequate servicing, or resulting from failure to comply with the recommendations contained in this manual. It is forbidden to use the product other than in accordance with its intended purpose!

NOTE: When a product-related "serious incident" occurs, which directly or indirectly led, could or could have led to any of the following events:

- a) death of a patient, user or other person or
- (b) temporary or permanent deterioration in the health of a patient, user or other person, or
- c) a serious threat to public health

the above "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTE: In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of a medical device, please consult a healthcare professional. **NOTE:** There may be a risk of overturning if the product is used improperly yourself.

NOTE: Please follow the sitting on the product/ standing up / moving instructions. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully. The maximum load of the rollator bag is 5 kg.

WARNINGS ASSEMBLY AND ASSEMBLY

After disassembling or reassembling the rollator, make sure that the seat is locked in both tubes and both grip tubes in the main support frame with bolts and nuts. Also make sure that the copper pins come out of the corresponding holes in the tubes of the wheels and that they are properly locked, and that both handles are at the same height and not moving.

WARNINGS INSTALLATION

Make sure the seat tube is securely fastened to the lock, all copper pins slide out of their respective holes, and the front and rear tubes are securely inserted into the main frame. BEFORE use, always check that the rollator and its components are securely locked in place.

UNFOLDING / FOLDING THE PADDLE

Place the rollator on a level surface, unfold the frame and wheels, then tighten the screws to align the front wheels and both tubes for the rear wheels (fig 1,2,3)

fig.1

fig.2

fig.3

ADJUSTING THE HANDLE TUBE

Press the bolt on the adjustment tube and slide it up or down to adjust the height. Make sure that the bolt is in the correct hole after adjusting.

FOLDING

Lift the handle of the seat to fold the rollator and then fold the wheels as shown in pictures 4 and 5

fig. 4


fig. 5

APPLICATION / LOCKING / UNLOCKING OF HAND BRAKES

CAUTION: When using and operating the rollator, and when folding it and adjusting the mechanisms, there may be a risk of entrapment and / or pinching of body parts of the user / accompanying person in openings / gaps between the components. These steps should be performed with extreme caution. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully.

APPLICATION OF HAND BRAKES

1. Pull the brake towards the handgrips.
2. Do one of the following actions:
 - A. STANDING - hold the brake handle up:
 - B. MOVING - Release the brake handle.

LOCKING THE HAND BRAKES

1. Push down on the bottom of the brake handle as shown in the picture until you hear a click (Fig. 7).

Fig. 6 (unlocked state)

Fig. 7 (locked state)

UNLOCKING THE HAND BRAKES

Pull the brake handles upwards until they release as shown in fig. 5.

1. Release the brake handle.
2. Make sure the brake cable is running along the side-line and the unused part is coiled in front of the rollator.

INTENDED USE: The product is intended to support the movement of people with reduced walking ability. The rollator is designed for one person.

CONTRAINDICATIONS: physical or mental limitations (e.g. visual impairment) that prevent safe handling with the rollator.

CHARACTERISTICS:

- Lightweight aluminium
- Ability to quickly unfold / fold the rollator
- Equipped with an extra bag
- Has a seat function

MAINTENANCE

1. Please ensure that all attachments are always secured.
3. Replace any damaged or destroyed parts immediately.

WHEEL MAINTENANCE

1. The wheels should be cleaned with warm water and a mild detergent. Pat dry with a dry cloth.
2. Lubricate periodically with 3-in-1 grease if the wheels begin to squeak.

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the medical device can be disposed of as normal household waste, except for electrical products - proceed in the manner appropriate for the disposal of electrical and electronic equipment.

Technical specifications:

Foldable, four-wheeled aluminium rollator.

Height: 79.5 - 92.5 cm

Seat height: 53.5 cm

Seat length: 46 cm

Seat width: 25 cm

Weight: 7.35 kg

Distance between the handles: 46.5 cm

Distance between the front wheels: 59 cm

Wheel diameter: front 8 "

Distance between the front and rear wheels: 51 cm

Maximum user weight: 136 kg

INCLUDED: ROLLATOR, OPERATING MANUAL

THIS MARK MEANS THE MAXIMUM WEIGHT OF THE USER

TARGET GROUP OF PATIENTS

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/size/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol.
7. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
8. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE), IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

Date of issue of the manual: October 10, 2023
v1-2 10/10/2023

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ENG

GENERAL NOTICE

Do not use the product until you have read and understood these instructions. If you do not understand the warnings, precautions, or instructions, contact your doctor or distributor to prevent damage to the product.

PURPOSE

The bag is a product that serves as an aid during transport.

It is forbidden to use the product for any purpose other than that for which it was intended.

CONTRAINDICATIONS

Physical or mental limitations (e.g., poor eyesight) that prevent safe handling of the product.

CHARACTERISTICS

- equipped with an extendable handle
- Castors for easy transport of the bag

The black transport bag from ANTAR can be easily opened by a long zipper, pack the folded rollator into it and close the bag. The bag has two handles for easy carrying and lifting, for example in the trunk of your car or as hand luggage on the plane.

This bag also has two wheels and on the other side an extendable n telescopic handle, so you can easily pull the bag behind you at the airport or train station.

The transport bag was specially developed for the travel rollator AT51006.

PARAMETER

Length: 66 cm

Width: 35 cm

CLEANING

The product should be cleaned with a damp cloth and a mild detergent. Do not dry clean, do not bleach, do not tumble dry.

STORAGE

The product should be stored in a dry, cool room, protected from moisture and direct sunlight.

HOW TO DISPOSE OF THE PRODUCT AT THE END OF ITS LIFE CYCLE

At the end of its service life, the medical device can be disposed of as normal household waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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Antar Medizin GmbH