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Introduction

Thank you for purchasing our product. Please read this user manual carefully before using it. Incorrect use may be harmful to health, and lead to injury or accidents. Thus, for safe and successful use, please read this manual carefully.

GENERAL WARNINGS

DO NOT assemble or use the product before reading and understanding this manual. If you do not understand the warnings, cautions or instructions, contact your healthcare professional, distributor or technician before attempting to assemble the product to avoid injury or damage to the product. Inspect all parts for damage caused in transport. If such damage is noticed, DO NOT use the product. Contact the manufacturer for more information.

NOTE: In the event of a device-related "serious incident" that has directly or indirectly led to, could have led to, or is likely to lead to any of the following:

- (a) death of the patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTE: In the event of pain, allergic reactions or other distressing symptoms, unclear to the user, associated with the use of the medical device, consult a healthcare professional.

NOTE: The manufacturer cannot be held responsible for damage caused by negligent maintenance, inadequate servicing or resulting from failure to follow the instructions in this manual.

ATTENTION: It is forbidden to use the product other than in accordance with its intended use!

CAUTION: There may be a risk of tipping over if the product is used incorrectly. Please observe the recommendations for getting on/off. It is forbidden to rock on the chair, as this may lead to the user tipping over and serious injury.

CAUTION: During use and handling of the product and when folding and adjusting the mechanisms, there may be a risk of entrapment and/or compression of parts of the user's/companion's body in openings/gaps between components. Perform these operations with particular care. Once the adjustment has been completed, stabilise the position by carefully tightening the nuts/bolts.

PURPOSE OF THE PRODUCT

The toilet chair is intended for disabled persons with musculoskeletal disorders and the elderly. It serves as a support during daily toileting.

CONTRAINDICATIONS

Physical or mental limitations that prevent safe handling of the product.

PRODUCT FEATURES:

- comfortable seat
- lightweight yet robust construction
- height adjustable

ASSEMBLING THE PRODUCT AT51009

Screw legs A and B to the seat, arranging them in an x-shape, using the five screws and washers provided, screwing them into the correct holes. Finally, adjust the seat height accordingly (6 degrees of adjustment in 2.5 cm increments) and check that the chair is stable and all screws are properly tightened.

ASSEMBLING THE AT51010

Screw the backrest tubes with 4 screws into the corresponding holes. Then fix legs A and B by arranging them in an x-shape and screw them to the seat. Then place the backrest tubes into the corresponding holes on the A and B frames until the locks lock (distinctive sound is heard). Finally, adjust the appropriate height of the seat (6 degrees of adjustment every 2.5 cm) and check that the chair is stable and all screws are properly tightened.

DISASSEMBLY OF THE PRODUCT

Perform all operations in the reverse order of the steps described in the "assembling the product" section.

STORAGE:

Store the product in a dry room at room temperature. Do not expose to direct sunlight. Do not store in areas with high humidity.

MAINTENANCE

1. Check all fixing screws before each use and tighten if necessary.
2. Promptly replace any damaged or deteriorated components.

CLEANING

Clean the product with a damp cloth and mild detergent. Do not use rough tools or strong detergents for cleaning. The product should then be wiped dry.

DESCRIPTION OF THE COMPONENTS AT51009

1. Plastic seat x1
2. Aluminium leg frames (A and B) x2
3. screws and washers to fix the seat x5

TECHNICAL SPECIFICATIONS AT51009

Overall width: 48 - 53 cm
 Height from floor to seat: 37 - 49.5 cm
 Overall length: 48 cm
 Seat width: 49 cm
 Seat depth: 29 cm
 Weight: 2.2 kg
 Maximum user weight: 120 kg
 Frame: aluminium
 Seat and backrest: PE

DESCRIPTION OF COMPONENTS AT51010

1. Plastic backrest x1
2. Plastic seat x1
3. Aluminium support frame connecting backrest and seat x1
4. Aluminium leg frames (A and B) x2
5. Bolts and washers for fixing the seat and back x5

TECHNICAL SPECIFICATIONS AT51010

Overall height: 69 - 81.5 cm

Overall width: 48 - 53 cm

Overall length: 48 cm

Seat width: 49 cm

Seat depth: 29 cm

Height from floor to seat: 37 - 49.5 cm

Weight: 3 kg

Maximum user weight: 120 kg

Frame: aluminium

Seat and backrest: PE

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER.

DELIVERY RANGE: chair (see description of components), user manual

TARGET GROUP OF PATIENTS

Persons suffering from diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this product is intended (see section on intended use of the product in this manual). The device can be purchased by the user independently or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/sizes/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

HOW TO DISPOSE OF THE DEVICE

The medical device can be disposed of as normal municipal waste.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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