AT51029, AT51027
INSTRUKCJA OBSŁUGI, CHODZIK
INSTRUCTION MANUAL, ROLLATOR
GEBRAUCHSANWEISUNG, ROLLATOR
NÁVOD K POUŽITÍ, CHODÍTKO
NÁVOD NA POUŽITIE, CHODÍTKO
MANUEL D'INSTRUCTION, DÉAMBULATEUR
GEBRUIKSAANWIJZING, ROLLATOR
MANUAL DE INSTRUCCIONES, ANDADOR
MANUALE DI ISTRUZIONI, DEAMBULATORE
BRUKSANVISNING, RULLATOR

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Thank you for purchasing our product. Before using the rollator, read the manual and instructions for assembling the product. Please follow the safety warnings.

GENERAL WARNINGS

To maintain safety while using the product, you **MUST** follow the guidelines below:

NOTE: DO NOT assemble and use the product before reading and understanding this manual. If you do not understand any warnings, cautions, or instructions, contact your healthcare professional or retailer before attempting to assemble the product to avoid injury or damage to the product. After each height adjustment, check the brake cable to avoid cable slack. Make sure that the knobs are locked in such a way that the main supporting structure will not move under the pressure. Failure to follow the above instructions may result in injury. Before use, check that ALL handle adjustments and that the wheels and moving parts of the product are in good condition. During use, ALL the wheels must **STAY** in contact with the ground at **ALL TIMES** to ensure the stability of the rollator. When the rollator is used in a stationary position, the hand brakes **MUST** be in the lock position.

NOTE: In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to, any of the following:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration of the patient's, user's or other person's health or
- (c) a serious threat to public health

this "serious incident" should be reported to the manufacturer and to the competent authority of the Member State where the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTE: In the event of pain, allergic reactions or other disturbing symptoms, unclear to the user, related to the use of the medical device, consult a healthcare professional.

NOTE: The manufacturer is not responsible for damage caused by neglected maintenance, improper servicing or failure to comply with the recommendations contained in this manual.

ATTENTION: It is forbidden to use the product in a way other than for its intended purpose!

ATTENTION: During the use and handling of the product, as well as during its assembly and adjustment of mechanisms, there may be a risk of entrapment and/or squeezing of the user's/assistant's body parts in the openings/gaps between the elements. These steps should be performed with particular care.

NOTE: If the product is used incorrectly, there may be a risk of tipping over. Please follow the getting on/getting off/moving instructions. After adjustment, stabilize the position by carefully tightening the nuts/bolts.

APPLICATION

A rollator is a device that facilitates the movement of patients with disabilities or who have difficulty moving. It can be used as an aid for patients after surgical procedures who do not show coordination disorders, patients with slight dysfunction of the lower limbs or their weakness. The rollator can also be used as an aid in moving geriatric patients.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

UNFOLDING

To unfold the rollator, slide the side frames outwards, then lock the horizontal expansion hinge. Then place the grip tubes in the correct holes and attach the brake cable to the frame. Finally, install the basket (assembly description in the instructions below).

FOLDING

All the steps described above should be performed in reverse order.

ASSEMBLY OF THE BASKET

Remove the tray from the basket before attaching the basket to the support. The basket should be placed on the hooks located at the height of the knobs. Then place the tray on top of the basket so that the cut ends are facing the handles.

NOTE: Before each use, check the operation of the brake system. Squeeze the brake handles to check that the brake blocks are in good contact with the wheels and are able to stop the rollator. Press the handles down to activate the handbrake and check that it works properly.

HANDLE HEIGHT ADJUSTMENT.

Handles should be installed in the openings of the side frames and adjusted to the correct height. A screw should be inserted from the inside of the side frame, and a washer should be put on the screw from the outside, and then the knob should be tightened. It is necessary to lock the frame with a knob. Screw the knob firmly to the support frame.

SAFETY WARNINGS ASSEMBLY AND ASSEMBLY

After unfolding or assembling the rollator, make sure that both handles are at the same height and do not move.

USING THE HANDBRAKES

1. Pull the brake towards the hand grips.
2. Do one of the following:
 - A. STOP - keep the brake handle up;
 - B. WHEN MOVING - release the brake handle

LOCKING THE HANDBRAKES

1. Push down the bottom of the brake handle as shown until you hear a click (Fig.6)
2. Refer to the UNLOCKING BRAKES section of this manual for more information


Fig. 5 (unlocked state)


Fig. 6 (locked state)

RELEASING THE HANDBRAKES

Pull the brake handles up until they release as shown in Fig. 5.

2. Release the brake handle.
3. Make sure the brake cable runs along the side tube and the unused part is rolled up at the front of the rollator.

TECHNICAL PARAMETERS:

AT51027:

Frame: steel
 Wheel size: 8"
 Weight: 6.5 kg
 Maximum width: 63 cm
 Maximum length: 65 cm
 Height: from 83 cm to 97 cm
 Width between handles: 26 cm
 The distance between the handles: 45 cm
 Tray and basket length: 28 cm
 Tray and basket width: 24 cm
 Basket height: 19 cm
 Bag height: about 35 cm
 Maximum basket load: 5 kg
 Maximum user weight: 120 kg

AT51029:

Frame: aluminum
 Weight: 6 kg
 Wheel size: 8"
 Maximum width: 63 cm
 Maximum length: 65 cm
 Height: from 83 cm to 97 cm
 The distance between the handles: 45 cm
 Tray and basket length: 28 cm
 Tray and basket width: 24 cm
 Basket height: 19 cm
 Bag height: about 35 cm
 Maximum basket load: 5 kg
 Maximum user weight: 120 kg

**THIS MARK INDICATES THE
MAXIMUM USER WEIGHT.**

MAINTENANCE

1. Ensure that all attachments are always secured.
2. Check the fixing screws of the handles. Tighten if necessary.
3. Promptly replace any damaged or damaged components.

WHEEL MAINTENANCE

1. Clean the wheels with warm water and a mild detergent. Dry with a dry cloth.
2. Lubricate periodically with 3-in-1 grease if wheels start to squeak.

TARGET GROUP OF PATIENTS

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

DELIVERY RANGE

Rollator, bag, basket, 2x knobs, 2x screws, 2x caps, instruction manual.

CLEANING

Clean the product with a damp cloth, do not use strong detergents or sharp objects.

STORAGE

The product should be stored in a dry, cool room, away from moisture and direct sunlight.

METHOD OF DISPOSAL OF THE PRODUCT AFTER WITHDRAWAL FROM USE

At the end of life, the medical device can be disposed of as normal household waste.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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