AT51040

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USER MANUAL ROLLATOR

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NÁVOD NA POUŽITIE ROLÁTOR

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INTRODUCTION

The rollator is manufactured from high quality materials and come with a tray, bag and locking brakes. The handles are adjustable to the correct height and the rollator folds easily for transport or storage.

APPLICATION/PURPOSE

Designed to assist a person with reduced mobility when moving and carrying small purchases. The user should have the necessary physical, visual and cognitive skills to operate the rollator safely. If in any doubt, consult your dealer. It is important that the user reads and understands the safety instructions in this manual. The product can be used repeatedly.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

-Lack of spinal stability: People with severe spinal stability problems or spinal injuries may not be suitable candidates for using a rehabilitation walker, as this can increase the risk of further injury.

-Uncontrolled seizures: In people with uncontrolled seizures or epilepsy, the use of a rehabilitation walker can be dangerous for both the patient and those around them.

-Poor muscle strength: In people with very limited muscle strength, using a rehabilitation walker may be difficult or ineffective. In such cases, a different approach to rehabilitation may be necessary, such as manual therapy or Muscle Strength Support Exercises.

-Balance problems: In patients with severe balance problems, the use of a walker can be dangerous, as it can increase the risk of falls and further injury.

-Rehabilitation walker mismatch: Improper fit of the rehabilitation walker to the patient can lead to additional injuries or worsening of the condition

USE


1. always make sure the brakes are in the locked position if the rollator is left unattended
2. check that all wheels are well assembled before use
3. check that the brakes function correctly
4. check that the handle height adjustment knobs are tightened before use
5. check that the rollator is fully unfolded and assembled
6. take care when carrying the folded rollator
7. keep your hands away from the folding mechanism to avoid the risk of pinching.

PRECAUTIONS

1. Always walk with the rollator in a smooth, controlled motion. Ask for assistance if necessary.
2. Do not hang objects on the brake levers or hand grips. This may affect the operation of the brake and/or the stability of the rollator.

LIST OF COMPONENTS

- 1) Tray
- 2) Handle
- 3) Two rear wheels
- 4) Main frame
- 5) Front wheels
- 6) Two front screws
- 7) Two T-bolts

SPECIFICATIONS

Overall width: 60 cm
 Overall length: 58 cm
 Tray height: 67 cm
 Tray size: 39.5*27.5 cm
 Adjustable handle height: 85-97 cm
 Weight: 6.2 kg
 Package size: 51*17*63 cm


Solid wheels: 6"
Max user weight: 108 kg

This symbol indicates the maximum weight of the user!

ASSEMBLY INSTRUCTIONS

- 1 Remove all parts.
2. Remove the main frame and open it as shown by the arrow. (1)
3. Insert the front wheels into the main frame. (2)
4. Insert the rear wheels into the main frame. (3)
5. Open the two 'on-off' buttons (4) and then insert the handle into the main frame and lock the 'on-off' button.
6. Place the tray on the main frame (5)

DISASSEMBLING

Perform all the above steps in reverse order.

ADJUSTABLE HANDLE

The handles of the rollator can be adjusted to suit different users. Simply open the two switches and set them to the desired height. (4) Please make sure that the two sides are at the same height.

BRAKE OPERATION ADJUSTMENT

When moving, the brake lever handle is in the middle position. Pull up on the brake handle lever to brake. Push down on the brake handle lever to lock the rear wheels. To release the lock, gently pull up on the brake handle lever so that it returns to the centre position.

Adjust the brakes as follows:

Turn the adjustment screw counterclockwise to increase the amount of pressure. Then lock the adjustment lock nut in a clockwise direction while holding the adjustment screw still. Brake adjustment reduced to light pressure may not be effective when braking. In this case, tighten the adjusting nut clockwise while holding the adjusting screw stationary.

MAINTENANCE

Over time, the brake cables can stretch and if no further adjustments can be made, the cables should be replaced by a service centre. Check that the wheels are secure each time before using the rollator. Also check that the brakes work properly and the handles and knobs are secure. If any part of the rollator appears loose or insecure, stop using it until an authorised dealer has carried out an inspection.

MAINTENANCE SCHEDULE:

Weekly

- Check brakes for malfunction
- Check all wheels and parts for loose bolts, nuts and knobs

Periodically

- Check hand grips for wear
- Check that frame bolts are tightened
- Check front and rear axle bolts are tightened
- Check material for wear or cracks

Annually

- Have the rollator serviced once a year for a full inspection, cleaning and lubrication of all moving parts.

CLEANING

The rollator can be cleaned with a damp cloth or sponge and warm soapy water. If necessary, gentle household antibacterial surface cleaners can be used. Take care not to leave any cleaning residue on the hand grips. It can make them slippery during use. Do not use scrapers or rough brushes as this may damage surface coatings. Do not use strong detergents for cleaning.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this device is intended (see the section on the purpose of the device of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer, should be taken into account.

WARNINGS!

- 1) Do not stand on this product.
- 2) Make sure both push handles are set at the same height.
- 3) Always keep a safe distance between your rollator and pedestrians or objects.
- 4) Do not attempt to replace the wheels of your rollator yourself. If they are faulty or need replacing contact your dealer.
- 5) Make sure the rollator is fully unfolded and assembled before using it again.

NOTES:

In the event of a device-related "serious incident" that has directly or indirectly led to, could have led to or may lead to any of the following:

- a) death of the patient, user or other person, or
- b) temporary or permanent deterioration of the health of a patient, user or other person, or
- c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

In the event of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

CAUTION: During use and handling of the device, and when assembling and adjusting the mechanisms, there may be a danger of the user's/companion's body parts becoming trapped and/or compressed in the openings/gaps between components. Perform these operations with particular care. Once adjustments have been completed, stabilise the position by tightening the nuts/bolts firmly.

ATTENTION: It is forbidden to use the product other than for its intended purpose.

NOTE: Before use, check that all components have been assembled correctly.

ATTENTION: The manufacturer is not liable for damage caused by negligent maintenance, inadequate servicing or as a result of failing to observe the recommendations contained in these instructions for use.

HOW TO DISPOSE OF THE PRODUCT AT END-OF-LIFE

After the device has been taken out of service, the medical device can be disposed of as normal household waste.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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