AT51056
INSTRUKCJA UŻYTKOWNIA, CHODZIK
INSTRUCTION MANUAL, ROLLATOR
GEBRAUCHSANWEISUNG, ROLLATOR
NÁVOD K POUŽITÍ, CHODÍTKO
NÁVOD NA POUŽITIE, CHODÍTKO
MANUEL D'INSTRUCTION, DÉAMBULATEUR
GEBRUIKSAANWIJZING, ROLLATOR
MANUAL DE INSTRUCCIONES, ANDADOR
MANUALE DI ISTRUZIONI, DEAMBULATORE
BRUKSANVISNING, RULLATOR

MD

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EN

Thank you for purchasing our product. Before starting use, please read this user manual carefully. Improper use may result in personal injury, damage, or an accident. Therefore, for safe and proper use of the product, please familiarize yourself thoroughly with this manual. **WARNING:** All parts must be checked for damage. Before first use, thoroughly inspect the technical condition of the product and ensure that all parts are complete and properly assembled (if applicable, pay particular attention to screws, nuts, rivets and other fastening and structural components). Make sure that the product is not defective or damaged. If any damage or missing parts are found, **DO NOT** use the product. For further information, contact the seller.

TARGET PATIENT GROUP

People suffering from diseases, dysfunctions, or injuries for the treatment, rehabilitation, or compensation of which this product is intended (see the section on the intended purpose of the product in this user manual). The product may be purchased independently by the user or on the recommendation of a physician, therapist, or other specialist. In both cases, the available sizes, required functions and variants of the product, indications and contraindications for use, as well as the information provided by the manufacturer, should be taken into account.

INTENDED PURPOSE

The product is intended to support mobility in people with limited walking ability. The rehabilitation rollator is designed for use by one person.

INDICATIONS

The product may be used in the following cases:

Balance training: The rehabilitation rollator enables patients to practice maintaining balance, which is essential in the rehabilitation process after injuries or neurological disorders.

Posture maintenance: Patients can use the rehabilitation rollator to train proper body posture, which is particularly important for people with postural defects or problems maintaining body stability.

Improvement of motor coordination: Exercises performed using the rehabilitation rollator may help improve motor coordination and movement control, which is important for people with neurological disorders.

Support in post-injury rehabilitation: The rehabilitation rollator may be used in the rehabilitation process after lower limb injuries, helping patients gradually regain mobility.

Increasing independence: For some patients, the rehabilitation rollator may serve as a tool to increase independence in daily life by enabling independent movement.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

Use of the rehabilitation rollator may be inappropriate for people with severe spinal instability or spinal injuries, as it may increase the risk of further injury.

In people with uncontrolled seizures or epilepsy, the use of a rehabilitation rollator may be dangerous for both the patient and people around them.

In people with very limited muscle strength, use of the rehabilitation rollator may be difficult or ineffective. In such cases, other rehabilitation approaches may be required, such as manual therapy or muscle-strengthening exercises.

In patients with severe balance disorders, use of the rollator may be unsafe as it may increase the risk of falls and further injuries.

LIST OF COMPONENTS

Main frame, handle with rear wheel (left), handle with rear wheel (right), backrest, shopping bag, front wheel, hex screw M6 – 3 pieces, five-point knob – 3 pieces, user manual.

PRODUCT FEATURES

Lightweight and easy to store or transport

Very good maneuverability and stability

Adjustable handle height for user adaptation

Equipped with a seat to allow resting

Includes a shopping bag for carrying small items

TECHNICAL SPECIFICATION

Frame: material and finish: aluminum, liquid coating

Seat: fabric

Backrest: Oxford fabric

Handles: PVC

Maximum user weight: 136 kg

Maximum bag load: 5 kg

Total weight: 5.2 kg

Length: 59.5 cm

Overall width: 71.5 cm

Overall depth: 59.5 cm

Folded dimensions (L x W x H): 69 x 25 x 79 cm

Seat width: 35 cm

Seat depth: 26 cm

Seat height: 55 cm

Distance between handles: 49 cm

Backrest height: 27 cm

Adjustable handle height (min–max): 79–101 cm

adjustment range: 25.4 cm

number of adjustment positions: 10

Swivel front wheels and rear wheels with brakes

Front and rear wheels: 8 inches
Handles equipped with reflectors
WARNING! Individual product dimensions may vary by ± 2 cm.


This symbol indicates the maximum user weight.

ASSEMBLY


Unfold the frame by pulling it apart to the sides.

Insert the handle with the rear wheel (left) into the main frame (arrow B).
Make sure that the side holes are aligned.
Insert the hex screw into the side hole of the main frame and push it all the way in (arrow C).
Fit the five-point knob onto the M6 hex screw and tighten it (arrow D).
Install the second handle in the same way (arrow E).
Assembly on both sides is identical.
Insert the front wheel into the main frame, aligning the holes (arrow F).
Insert the hex screw and push it fully in (arrow G).
Tighten the five-point knob.

Attach the backrest to the handles (arrows H/I).
Attach the shopping bag to the frame (arrows J/K).

FOLDING THE FRAME

Pull the strap under the seat to fold the frame. To unfold the frame, push the armrests outward.

HANDLES

The handles are height adjustable.

USING THE BRAKES

To stop the rollator, pull the brake levers upward toward the handles.

To remain stationary, keep the brake levers pressed down.

To start moving, release the brake levers – they will automatically return to their original position and the brakes will be released.

Always engage the brakes when you do not intend to move the device.

Never leave the device unattended without the brakes engaged when someone is seated on it.

Never leave the device on slopes without the brakes engaged.

When sitting on the rollator, always sit upright with your forearms or hands lightly resting on the handles.

Never lean backward while sitting on the rollator.

Locking / Unlocking the Hand Brakes

To lock the brakes, press the brake levers downward until you hear a “click.”

To unlock the brakes, pull the brake levers upward until you feel them release.

WARNING

In the event of a “serious incident” related to the product that directly or indirectly has led, could have led, or may lead to any of the following:

- a) death of a patient, user, or another person,
- b) temporary or permanent deterioration of the health of a patient, user, or another person, or
- c) a serious public health threat,

the incident must be reported to the manufacturer and the competent authority of the EU Member State in which the user or patient resides. If pain, allergic reactions, or other worrying or unclear symptoms occur during use of the medical device, consult a healthcare professional. It is prohibited to use the product in any manner other than its intended purpose.

CLEANING: Clean the device regularly and keep it clean. Light dirt and dust can be removed with a damp cloth. Heavier dirt can be removed using a mild, non-corrosive detergent and warm water. Do not use abrasive cleaners or chlorine-based cleaning agents, as they may cause damage or discoloration. Dry all parts with a soft cloth.

Wheels: Clean the wheels with warm water and a mild detergent. Dry with a clean cloth. Lubricate periodically with a multi-purpose oil or if the wheels start to squeak.

MAINTENANCE

To ensure safe and proper use of the product, regularly check individual components according to the following schedule: Before each use, check the technical condition of the product for damage and missing parts.

Weekly inspection: Check that screws are properly tightened. Inspect brake cables for looseness or damage. Check bolts and nuts and tighten if necessary.

Monthly inspection: Check handle grips for stability and damage. Ensure the brake levers function correctly. Check the bases (legs) for stability and looseness.

Quarterly inspection (every four months): Inspect the wheels to ensure they rotate freely and are undamaged. Check the bag for damage and usability. Inspect the frame (frame and supports) for cracks or deformation.

TRANSPORT: The rollator is easy to fold, making it convenient for transport. When transporting, it is recommended to fold it as described above. The rollator should be transported in a manner that ensures stability and prevents sudden or dangerous movements.

REUSE: The rollator is suitable for reuse. Before being passed on to another user, it must be inspected by a qualified dealer. The following points must be checked: Visual condition of the Frame Secure fastening of bolts and screws Condition of plastic components Condition of wheels and ball bearings Effectiveness and condition of the brakes If any component is damaged or broken, it must be replaced. Wheel replacement is recommended whenever any changes are noticed. All loose screws and nuts must be tightened. The product must be thoroughly cleaned and disinfected before reuse.

WARNING: Reuse of the rollator is prohibited if the frame is bent or broken. The user manual must accompany the rollator whenever the user changes. **STORAGE:** Store the product in a dry place and do not expose it to direct sunlight. **DISPOSAL:** After the product is withdrawn from use, it may be disposed of as regular household waste. Please follow local waste segregation regulations and protect the environment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawiślńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawiślńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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