AT51113
CHODZIK INSTRUKCJA OBSŁUGI
ROLLATOR USER MANUAL
GEBRAUCHSANWEISUNG FÜR GEHHILFEN

PL

Wstęp

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ENG

Introduction

Thank you for purchasing our product. Before using it, please read this user manual carefully. Incorrect use may lead to damage to health, injury or accidents. Thus, in order to use the product safely and successfully, please read this manual carefully.

FEATURES

The AT51113 rollator is made of high-quality, durable aluminium. The folding frame makes transport easier. The rollator is equipped with swivelling front wheels and non-twisting rear wheels with a brake and a parking clip. The product is enriched with a seat and backrest, making it possible to rest in a sitting position. The rollator also has a convenient bag for small shopping.

PURPOSE/ INDICATIONS

The rollator is intended for people who need support in moving around. It is suitable for people who want to take an active part in everyday life.

CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent safe handling.

WARNINGS

Check the condition of the product before carrying out any operation. If the product malfunctions or looks abnormal, immediately contact the dealer from whom you purchased the unit.

GENERAL PRECAUTIONS FOR USE

Always ensure that the handles are at the correct height and that they are properly secured before using. To avoid the risk of falling, the rollator must not be used on uneven or sloping surfaces. Always check the condition of the wheels, brakes and overall stability of the frame before use.

CAUTION: When using and operating the rollator and when folding and adjusting the mechanisms, there may be a risk of the user/assisted person's body parts becoming trapped and/or pinched in the openings/gaps between components. Perform these operations with particular care. Once the adjustment is complete, stabilise the position by tightening the nuts/bolts carefully.

CAUTION:

In the event of pain, allergic reactions or other distressing symptoms unclear to the user related to the use of the medical device, a healthcare professional should be consulted.

CAUTION:

In the event of a device-related "serious incident" that has directly or indirectly led to, could have led to, or is likely to lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

ATTENTION: The manufacturer shall not be held liable for damage caused by negligent maintenance, inadequate servicing or resulting from failure to observe the recommendations contained in this instruction manual.

ATTENTION: It is forbidden to use the product other than in accordance with its intended use!

ATTENTION: There may be a risk of tipping over if the product is used incorrectly. Please observe the recommendations for getting on/off/moving. After adjustment, stabilise the position by tightening the nuts/bolts carefully.

SCOPE OF DELIVERY

Walker, side frames (x2), armrests with brakes (x2), cane support, bag, user manual.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

PREPARATION FOR USE

The selection of the type of assistive device and its appropriate configuration must be made by qualified medical and technical personnel based on the medical conditions and the conditions of use of the device by the user.

ASSEMBLING


1. Unfold the frame by lowering the transverse seat tubes, then press down on the sides of the seat until the tubes lock and a distinctive sound is heard. Make sure the walker is correctly unfolded.
- 2 Fit the two side frames with the copper locks and adjust to the correct height. (1/2)
3. the armrests, including the brakes, should be inserted into the holes in the side frames and tightened using the knobs. (3)
4. Finally, fit the cane holder and attach the bag. (5)

ASSEMBLING

Follow the steps described above in reverse order.

USING THE BRAKES

- 1 To release or stop the rollator, pull the brake lever. Release the brake lever to move again (Fig. 6).
2. To activate the parking brake, push the brake lever downwards until it locks (fig. 7).

SPECIFICATIONS

total width 63 cm
width between handles 50 cm
total depth 71 cm
width of seat 46 cm
depth of seat 25 cm
height of seat 54 cm
height adjustment range of handles 97 - 127 cm
castors 8"
weight 8.4 kg
max. user weight 136 kg

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

STORAGE AND DISPOSAL

Store the baby rollator away from sources of light, heat and moisture. Avoid exposing the rollator to direct sunlight. Store in a dry and cool place. Dispose of the rollator in accordance with current environmental and recycling regulations.

MAINTENANCE AND CLEANING

After each use, remove accumulated dirt using mild cleaning agents, preferably soapy water and then wipe off with a dry cloth. Do not use sharp, abrasive objects. Periodically check the function of the brakes, the condition of the seals of the attachment parts and the condition of the welds.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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