AT51115

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Thank you for purchasing our product. Please carefully read the user manual. Do not use the product before familiarizing yourself with this manual and understanding its contents. If you do not understand the warnings, notes, and recommendations, contact a healthcare professional or the seller to avoid product damage or personal injury.

WARNING: Check all parts for damage incurred during transportation. If such damage is detected, DO NOT use the product. More information is available from the seller.

TARGET PATIENT GROUP

Individuals dealing with diseases, dysfunctions, or injuries for the treatment/rehabilitation or compensation of which this product is intended (see the section regarding the product's purpose in this user manual). The product can be purchased independently by the user or based on the recommendations of a doctor, therapist, or other specialist. Regardless of whether the product is purchased independently or recommended by a doctor/therapist/other specialist, consideration must be given to the available sizes/necessary functions and product variants, indications, and contraindications for use, as well as the information provided by the manufacturer.

BEFORE USE, READ THE ENTIRE MANUAL CAREFULLY

Walking canes feature single-height adjustment by setting the distance from the handle to the support point with a foot tip (rubber cap). Individual products are adjusted to the user's height by telescopically sliding the tubes to the desired height and then locking them in place.

The product's height is considered correct if, during use, the user maintains an upright posture and the cane operates vertically, i.e., parallel to the body's axis. Proper height adjustment and correct usage of our products significantly extend the cane's lifespan and provide the user with a sense of safety and satisfaction.

LIST OF COMPONENTS

Cane, user manual

WARNING:

In the event of a "serious incident" related to the product that directly or indirectly led, could have led, or may lead to any of the following events:

- death of the patient, user, or another person, or
- temporary or permanent deterioration of the health of the patient, user, or another person, or
- a serious threat to public health,

this "serious incident" must be reported to the manufacturer and the competent authority of the Member State where the user or patient resides. In Poland, the relevant authority is the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

WARNING:

If pain, allergic reactions, or other concerning, unclear symptoms associated with the use of the medical device occur, consult a healthcare professional.

WARNING: The product must not be used in any way other than its intended purpose.

FEATURES AND TECHNICAL SPECIFICATIONS AT51115(5)

Universal (can be used as a right or left version).

Made of aluminum.

Polypropylene handle.

Incremental height adjustment every 2 cm (maximum height 90 cm, minimum height 68 cm).

Maximum load: 120 kg.

Rubber tip: flexible, black (non-marking for floors and parquet).

Usage period: 5 years from the production date.

This symbol indicates the maximum user weight.

FEATURES AND TECHNICAL SPECIFICATIONS AT51115(1)/AT51115(3)/AT51115(3)

Universal (can be used as a right or left version).

Made of aluminum.

Plastic handle.

Incremental height adjustment every 2 cm (maximum height 92 cm, minimum height 70 cm).

Maximum load: 100 kg.

Rubber tip: flexible, black (non-marking for floors and parquet).

Usage period: 5 years from the production date.

This symbol indicates the maximum user weight.

APPLICATION:

For individuals who need support while walking, after lower limb injuries, to maintain balance, and to transfer body weight in traumatic and orthopedic conditions. Also for elderly people as an aid in maintaining balance while walking.

CONTRAINDICATIONS:

Coordination and balance disorders, upper limb injuries, and dysfunctions that prevent reliance on the cane.

WARNINGS

The cane is intended for use indoors and outdoors on even surfaces.

Use extreme caution when using the cane on wet, icy, greasy, or heavily sanded surfaces, as this increases the risk of falling.

Periodically inspect the cane visually, including checking the height adjustment mechanism and the condition of the tip. If any part of the cane is damaged, stop using the product immediately, as it poses a fall risk.

Protect from moisture!

Protect from direct sunlight!

Do not exceed the maximum allowable load!

Regularly inspect the product during and before use to ensure there are no scratches, cracks, or other damages.

USER INSTRUCTIONS

Loosen (turn counterclockwise, up to two full turns) the locking sleeve nut.

Press the height adjustment button with your thumb.

Slide the adjustable tube along the axis to the desired height, ensuring the holes in the adjustable tube align perfectly with the body holes.

Confirm that the adjustment button clicks into the correct hole (audible click).

Tighten the locking sleeve nut.


Note: Do not over-loosen the adjustable tube, as it may become improperly locked.

WARNING!


The handle of the cane should be at hip level (precisely at the greater trochanter of the femur), and the hand should be slightly bent at the elbow (no more than approximately 30 degrees). The height of the cane is correct if it remains vertical and parallel to the body axis during use.

The cane should be held in the hand opposite the injured or weaker leg. Move the cane at the same time the injured/weaker leg steps forward, while the body is supported by the healthy leg. The movement of the cane and the injured/weaker leg should be simultaneous. When the injured/weaker leg bears weight, part of this support should be transferred to the cane resting on the ground.

ADJUSTING CANE HEIGHT TO BODY HEIGHT:

CORRECT POSITIONING OF THE INDEX FINGER:

MAINTENANCE AND CLEANING

Clean the cane shaft and handle with a damp cloth using mild detergents.

Regularly remove dirt from the moving parts of the cane (height adjustment mechanism).

Do not use abrasive cleaners!

Do not clean the cane with pressure or steam cleaners!

WARNING: Do not use cleaning agents based on gasoline, solvents, acids, chlorine, etc.!

STORAGE

Store the product in a dry place, away from direct sunlight. Storing the product in damp environments (bathrooms, laundries) or outdoors can deteriorate the product's functional and aesthetic parameters.

DISPOSAL METHOD AFTER DECOMMISSIONING

After the product's life cycle, the medical device can be disposed of as regular municipal waste. Remember proper waste segregation to protect the environment.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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