AT51116
LASKA INSTRUKCJA OBSŁUGI
WALKING STICK OPERATING INSTRUCTIONS
GEHSTOCK GEBRAUCHSANWEISUNG
NÁVOD K POUŽITÍ
NÁVOD NA POUŽITIE
MANUEL D'INSTRUCTIONS
HANDLEIDING
MANUAL DE INSTRUCCIONES
MANUALE DI ISTRUZIONI
BRUKSANVISNING

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ENG

GENERAL WARNING

Do not use the product before reading and understanding this manual. If you do not understand any warnings, cautions, or instructions, contact your healthcare professional or retailer to avoid damaging the product.

NOTE: Inspect all parts for shipping damage. If you notice such damage, DO NOT use the product. More information from the manufacturer.

ATTENTION: it is forbidden to use the product in a way other than for its intended purpose

ATTENTION:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to, any of the following:

- a) death of a patient, user or other person
- b) temporary or permanent deterioration of the patient's, user's or other person's health or
- (c) a serious threat to public health

this "serious incident" should be reported to the manufacturer and to the competent authority of the Member State where the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

ATTENTION:

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of the medical device, consult a healthcare professional.

APPLICATION / INDICATIONS

The cane is intended for people who need support when walking, after lower limb injuries, to maintain balance and transfer body weight in traumatic and orthopedic diseases. Recommended for the elderly as an aid in maintaining balance while walking.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment, balance and coordination disorders, upper limb dysfunctions) that prevent safe handling of the product.

PRODUCT DESCRIPTION AND FEATURES

- lightweight, made of aluminum
- universal, can be used in the right and left version
- maximum load 120 kg
- non-slip overlay at the end
- folded into 4 elements
- adjustable height: 80-90 cm (5 degrees of adjustment every 2.5 cm)
- weight: 0.34 kg

USE

The product is packaged in a folded position. After removing from the package, unfold the cane, set it to the desired length.

CLEANING/MAINTENANCE

The product should be cleaned with a damp cloth with a mild detergent and then wiped with a dry cloth. Do not use abrasive objects, solvents or strong detergents.

Remove impurities from the moving parts of the cane on an ongoing basis (ratchet mechanism of the locks).

STORAGE

The product should be stored in a dry, cool room, away from moisture and direct sunlight.

TARGET GROUP OF PATIENTS

The health care professional points out the use of the device for adults and children at his own responsibility, taking into account the available variants/sizes/necessary functions/size and indications, taking into account the information provided by the manufacturer.

HOW TO DISPOSAL OF THE PRODUCT

After the device is decommissioned, the medical device can be disposed of as ordinary municipal waste

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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