AT51125, AT51127
INSTRUKCJA OBSŁUGI ROTOR REHABILITACYJNY, ELEKTRYCZNY
USER MANUAL PEDAL EXERCISER, ELECTRIC
GEBRAUCHSANWEISUNG
AT51125 ELEKTRISCHER REHABILITATIONSROTOR
MANUEL DE L'UTILISATEUR EXERCICE A PEDALE, ELECTRIQUE
GEBRUIKERSHANDLEIDING PEDAAL TRAINER, ELEKTRISCH
NÁVOD K OBSLUZE
REHABILITAČNÍ ROTOR, ELEKTRICKÝ
NÁVOD NA POUŽITIE
REHABILITAČNÝ ROTOR, ELEKTRICKÝ
MANUAL DE USUARIO
ROTOR DE REHABILITACIÓN, ELÉCTRICO
MANUALE D'USO
ROTOR RIABILITATIVO, ELETTRICO
BRUKSANVISNING
ELEKTRISK REHABILITERINGSROTOR

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EN

GENERAL WARNINGS

Failure to observe the following safety precautions may result in injuries that are not conducive to recovery or other injuries.

Warning:

- 1. Do not stand on the pedals.
- 2. Sit on a chair, sofa or bed when exercising.
- 3. Before using the product, please ensure that all bolts, screws, knobs, frames and brackets are fastened.
- 4. Patients, persons with disabilities and the elderly must use it accompanied by another person.
- 5. People with heart disease or who are prohibited from strenuous exercise should follow their doctor's advice and exercise accordingly.
- 6. This product should not be used by people with untreated injuries.
- 7. Children and people with disabilities must use the equipment under surveillance of guardians/carers.
- 8. This product is specially designed for family use.
- 9. In order to avoid injuries when transporting the equipment, please pay particular attention to safety issues.
- 10. Before using the product, perform appropriate stretching exercises to prevent injury during exercise.

CAUTION: Inspect all parts for any potential damage caused by transportation. DO NOT use the product in case of detecting any such damage. Contact the manufacturer for additional information.

CAUTION: Any serious incident involving the product should be reported to the manufacturer and the competent body of the Member State where the users or patient reside.

CAUTION: There may be a risk of entrapment and/or compression of parts of the body of the user/accompanying person in the openings between the components during use and operation of the product, and when assembling and adjusting the mechanisms. These operations must be carried out with particular care. After adjustment, stabilise the position by tightening the nuts/screws securely.

CAUTION: It is prohibited to use this product in violation of its intended purpose.

SAFETY MEASURES

Please read the following instructions carefully before using the product.

Consult your doctor before using the rotor or any other exercise device. This rule shall be observed in particular in the following cases:

- 1. Cancer.
- 2. Recovery after a recent surgical procedure.
- 3. Documented heart problems.

4. Issues with the hips.
5. Thrombosis.
6. Broken bones.
7. Epilepsy.
8. DO NOT USE during pregnancy.

CAUTION:

Periodically check that all fastening elements are properly tightened and properly connected. Every three months check and tighten all parts of the device.

If you have any doubts regarding any parts of the machine, please contact your dealer.

We recommend replacing worn components with original parts. The use of other spare parts may cause injury or affect the operation of the machine.

GENERAL INSTRUCTIONS

Certain precautions, including the ones listed below, should always be observed when using electrical appliances.

This includes important information on safety, use and maintenance.

USE

Rehabilitation rotors is a compact exercise bicycle which can be easily stored if needed. The rotor is used primarily for rehabilitation of people with limb paresis caused, for instance, by an accident or a stroke. When used regularly, it upregulates blood pressure and increases muscle strength. It is a great method for moderate cardiovascular exercises.

Connect the power supply plug to a DC power source and press the red button under the remote control to position I. When the button (1) is pressed, the unit will run at speed level 1 and time will be set to 5 min. The display is going to show cyclical random data: speed, time, distance, calories (scanning).

Pressing the (2) Time button allows you to set time ranging from 1 to 15 minutes. For faster adjustment, hold down the button (2) until the desired time is displayed.

The device switches off automatically after the set time has elapsed.

Press button (4) Speed+ to increase the speed. Press button (3) Speed- to reduce the speed. For faster adjustment, hold down button (4) or (3) until you reach the desired speed level.

Use the (5) Reverse button to change the direction pedals rotation.

Press (6) MODE at once to display speed, time, distance and calories. Press the button several times until a red light appears under the letters Scan in order to display random parameters. Press (1) to switch off. When finished, press the red button below the remote-control bracket. Setting it to 0 switches the unit off completely.

INSTRUCTION MANUAL

It is recommended to start with short exercise sessions of 3 to 6 minutes, at low speed, twice a day. Gradually increase the duration of sessions to a maximum of 15 minutes per session, twice a day. Once the user is comfortable with the exercise, try higher speeds. However, we recommend prior consultation with your doctor.

ASSEMBLY INSTRUCTIONS

Remove the unit from the box and lay it on the floor, as shown in Figure 1.

Assistance from another person is required for assembly work

Scope of delivery:

(A) Main body; (C-L) Left pedal (C-R) Right pedal; (B) Stabiliser bar; (B1) Rear stabiliser bar; (D) M-4x25 bolt; (E) M-4 lock washer; Pedal combination spanner.

1.INSTALLING STABILISERS

Take the stabiliser (B) marked (A) and place it on the front of the unit marked (A) as shown in Fig. 2, then insert the M-4x25 screws (D) together with the lock washers (E).

Next, take the rear stabiliser (B1) marked with the letter (B) and follow the same assembly procedure as for the front stabiliser, Fig.3,

2.INSTALLING PEDALS

Right and left mean the position the user takes when performing the exercises.

Insert the right pedal (C-R), marked with (R), onto the right crank, also marked with (R). Next, Fig.4. Install the left pedal (C-L), marked with (L), on the left crank, also marked with (L), Fig.4.

3.CONNECTION TO THE MAINS POWER SUPPLY

Connect the device to an earthed plug (220 V). Set the switch (L) in position [1] Fig.5.


Fig. 1


Fig.2

Fig.3


Fig. 4

Fig.5

AT51125

Display indicator: Speed, Time, Distance (distance in km), Calories, Scan (shows random samples of data, such as speed, time, distance, calories).

- 1 - Power supply
- 2 - (Time) Button for time selection (1-15 mins).
- 3 - Speed - button for reducing speed (12 speed levels).
- 4 - Speed - button for increasing speed (12 speed levels).
- 5 - (Reverse) Forward/reverse - button for changing the direction of pedal movement.
- 6 - (MODE) Button for selecting display mode of speed, time, distance and calories


REMOTE CONTROL

Scanning Speed Time Distance Calories

AT51127

7. on/off button
8. direction change
9. speed setting
10. time setting
- Change mode (active/passive) 12.
12. speed
- Distance
14. calories
15. training time

Weight: 5kg
Width: 30cm
Length: 50 cm
Height: 28 cm
Maximum user weight: 120kg

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

CAUTION! In order to reduce the risk of electric shock:

Unplug the unit from the power outlet immediately after use and before cleaning. Please carefully read the instructions contained in this document. It contains important information on safety, use and maintenance of the treadmill.

1. Make sure that the supply voltage is suitable for the machine (220 V), Fig.5 before connecting the machine. The supply

voltage must match the one indicated on the rating plate and the holes in the wall socket must match the pin configuration of the appliance plug. If this is not possible, consult your dealer.

2. This device has been designed exclusively for domestic use.
3. Do not leave the unit unattended while it is plugged in. Unplug the unit from the electrical outlet after each exercise session.
4. Place the unit on a flat surface and leave at least 0.5 metres of free area behind it.
5. Parents and/or carers of children must always consider their curious nature, which may often lead to dangerous situations and behaviours causing accidents. Therefore, they must always be supervised. The device must not be used as a toy under any circumstances.
6. Keep hands away from moving parts.
7. Wear appropriate clothing and footwear. Make sure all laces/strings are properly tied.
8. It is not advisable to store the unit in a damp place for a long time as it will be prone to rusting.
9. Do not use the unit if the power cord is damaged or worn.
10. Keep the power cord away from hot surfaces.
11. This machine is intended for outdoor use.
12. If you wish to disconnect the unit from the power supply, first turn the power off using the button and then remove the plug from the wall socket.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

CLEANING/MAINTENANCE

Turn off the unit and unplug the power cord before cleaning! Clean with a slightly damp cloth and soap.

In order to maintain the performance of the product, please pay particular attention to cleaning all the components of the product, and to store it in proper conditions. Storage in rooms with very high humidity (baths, laundries, etc.) or in the open air, combined with poor maintenance, will result in a rapid loss of performance and aesthetics.

CAUTION: Do not expose the product to direct sunlight as certain elements may discolour.

CAUTION: The manufacturer shall not be held responsible for any damage caused by negligent maintenance, improper servicing or being the result of a failure to observe the recommendations contained in this instruction manual.

WARNING:

In the event of a device-related “serious incident”, which directly or indirectly led to, could have led to, or is likely to lead to any of the following events:

- a) death of a patient, user or another person, or
- b) temporary or permanent deterioration of the health of a patient, user or another person, or
- c) a serious risk to public health, such a “serious incident” must be notified to the manufacturer and to the competent authority of the Member State, in which the user or patient resides. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

WARNING:

In the event of pain, allergic reactions or other alarming, unclear symptoms related to the use of the medical device, please consult a healthcare professional.

NOTE: It is forbidden to use the product contrary to its intended use!

NOTE: Inspect all parts for damage caused in transit. If any such damage is found, DO NOT use the product. More information from the manufacturer.

NOTE: Before use, check that all components are properly assembled.

TARGET GROUP OF PATIENTS

The healthcare professional indicates the use of the product to adults and children at their own risk, taking into account the available variants / sizes / necessary functions / size and indications, taking into account the information provided by the manufacturer.

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the electrical products - proceed in the manner appropriate for the disposal of electrical and electronic equipment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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