AT53029

**INSTRUKCJA OBSŁUGI KOŁNIERZ SZYJNY MIĘKKI
USER MANUAL SOFT CERVICAL COLLAR
BENUTZERHANDBUCH WEICHE HALSKRAUSE
NÁVOD K OBSLUZE MĚKKÝ KRČNÍ LÍMEC
NÁVOD NA POUŽITIE MÄKKÝ KRČNÝ GOLIER
MANUEL D'UTILISATION COLLIER CERVICAL SOUPLE
GEBRUIKERSHANDLEIDING ZACHTE CERVISCHE KRAAG
MANUAL DEL USUARIO COLLARÍN CERVICAL SUAVE
MANUALE D'USO COLLO CERVICALE MORBIDO
ANVÄNDARHANDBOK MJUK CERVICAL KRAGE**

PL

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ENG

INTENDED USE / INDICATIONS

- diseases of the cervical spine, which lead to various deformations of the body posture,
- incorrect head positioning, in which the use of the collar will contribute to its gradual correction,
- degenerative changes in the cervical spine, i.e. discopathy, scoliosis,
- injuries as a result of accidents, contusions, injuries
- to relieve the joint parts, after operations, to stiffen the spine and thus prevent its breach
- as an aid during rehabilitation, it is also very useful when moving or transporting the patient

The period of wearing a brace is usually determined by the attending physician and may vary depending on the type of injury or the severity of the disease.

CONTRAINdications

physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

WARNING

Do not use the product after applying warming ointments or gels.

CHARACTERISTICS

Flesh color. Fastened at the back with Velcro.

WASHING INSTRUCTIONS

Hand wash in lukewarm soapy water. Do not use bleach and strong detergents. Do not iron or dry clean.

SIZES

In order to find the correct size, measure the circumference of the neck.

Size	S			M			L			XL		
Height	8	10	12	8	10	12	8	10	12	8	10	12
Neck circumference	<38 cm			38-42 cm			44-48 cm			49-53 cm		

CAUTION

When a product-related "serious incident" occurs, which directly or indirectly led, could or could have led to any of the following events:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person, or
- c) a serious threat to public health

this "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of a medical device, consult a healthcare professional.

ATTENTION

it is forbidden to use the product other than in accordance with its intended purpose

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the medical device can be disposed of as normal municipal waste, with the exception of electrical products - proceed in the manner appropriate for the disposal of electrical and electronic equipment.

TARGET GROUP OF PATIENTS

The healthcare professional indicates the use of the product to adults and children at their own risk, taking into account the available variants / sizes / necessary functions / size and indications, taking into account the information provided by the manufacturer.

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WARRANTY CARD

Model:
Serial number:
Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

Date of issue of the manual: September 16, 2022
v1- 26/09/2022

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