AT53031

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
USE OF THE PRODUCT

1. Choose the correct size for the brace by measuring the circumference of the wrist.
2. Adjust the Velcro strap of the brace to make it easier to insert the hand inside the brace.
3. Slip the brace over the hand.
4. Adjust the Velcro fasteners to ensure the wrist is adequately stabilised.
5. The brace has a removable reinforcement splint which can be used if required.

SIZE CHART (measure wrist circumference to find the correct size)

- S (13-15 cm)
- M (15-18 cm)
- L (18-20 cm)
- XL (20-23 cm)

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling.

INDICATIONS:

- strengthening and stabilising a weakened wrist;
- carpal tunnel syndrome;
- dislocations of the wrist joint or overuse syndrome;
- after wrist surgery;

FEATURES:

- wrist orthosis with a universal cut for the left and right hand
- an aluminium profiled insert and an extended strap system ensure complete support of the wrist, providing freedom of movement of the fingers
- made of breathable material warms the wrist and improves blood circulation

CAUTION:

In the event of a device-related 'serious incident' that has directly or indirectly led to, could have led to or is likely to lead to any of the following:

- (a) death of the patient, user or other person or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTE: In the event of pain, allergic reactions or other distressing symptoms, unclear to the user, associated with the use of the medical device, consult a healthcare professional.

NOTE: It is forbidden to use the product other than in accordance with its intended use

ATTENTION: The manufacturer will not be held responsible for any damage caused by negligent cleaning or as a result of not following the instructions in this manual.

WASHING INSTRUCTIONS

Remove metal parts, fasten all Velcro fasteners before washing. Hand wash in lukewarm water using mild soap. Do not dry clean or tumble dry. Do not use bleach. Do not expose to direct sunlight.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.


TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/size/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

DISPOSAL OF THE DEVICE AT END-OF-LIFE

At the end of life, the medical device can be disposed of as normal household waste.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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