## AT53037 AT53049

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**PURPOSE**

The ankle brace is designed for people whose ankle joint is weakened. Thanks to the cross straps, it is possible to adjust the range of flexion and extension movement of the foot. The brace is designed to prevent injuries to the joint in the post-traumatic period and to help stabilise the ankle joint. The brace is equipped with elastic inserts which help stabilise the ankle joint during lateral and torsional movements.

**APPLICATIONS:**

- Serves as a preventative measure to protect against injury and trauma.
- The brace with a system of reinforcements is suitable during rehabilitation, after sprains of the joint or after twisting the ankle, or immediately after removing the plaster cast.
- For active athletes when resuming training after a break caused by an injury or as a preventative measure to protect the joint from overload injuries.

**CONTRAINDICATIONS**

Physical or mental limitations (e.g. visual impairment) that prevent safe handling.

**RULES OF USE**

1. It is recommended that the indications and rules of using the orthosis are advised by a physician.
2. The method of putting on the orthosis should be demonstrated by qualified personnel, medically authorized to do so, in the presence of the person who takes care of the user (if applicable).
3. The brace is a Class I medical device intended for one person, therefore it cannot be rented to other users, even for a short period of time.
4. Use the product and its components only as intended.
5. The orthosis should be kept clean in accordance with the instructions given by the Manufacturer in this manual and on the label placed on the product. Improper cleaning and maintenance of the product (e.g. too rare washing, insufficiently thorough rinsing) may lead to skin disorders.
6. The skin in the place where the orthosis adheres to the body should be clean - while wearing the brace, do not apply therapeutic compresses made of ointments and creams, do not apply them just before putting the brace on.
7. Continuous use of the orthosis for a long period of time may result in muscle atrophy, therefore it is recommended to systematically perform exercises selected by qualified medical personnel under the supervision of the physician.
8. The wearing time of the orthosis should be determined by a specialist.

**FEATURES:**

- a unique system of reinforcements to increase the stability of the joint
- flexible stays that help stabilise the ankle joint during lateral and torsional movements
- fixing straps for perfect support of the joint
- light, air-permeable materials and low profile of the orthosis allow its use with any type of footwear.

**CLEANING AND MAINTENANCE**

1. The manufacturer recommends washing the product periodically
2. Close all velcro straps of the orthosis before washing
3. Remove metal parts before washing
4. Hand wash in lukewarm water with soap or neutral, anti-allergic detergent only, up to 40°
5. Rinse with lukewarm water. Do not use fabric softeners for the last rinse.
6. Do not wring out - remove excess water by hand pressure or in a towel
7. Do not bleach
8. Do not tumble dry
9. Dry the orthosis flat in a shaded and airy place - do not expose the orthosis to direct sunlight; do not use external heat sources for drying indoors
10. Do not iron
11. Do not dry clean
12. Do not use professional water cleaning (service laundries)

**STORAGE**

Store the washed and completely dried product folded, in a ventilated place without direct sunlight. Do not store in humid places or near artificial sources of heat. Store at room temperature.

### **HOW TO DISPOSE OF THE PRODUCT AT END-OF-LIFE**

Once the device has been withdrawn from service, the medical device can be disposed of as normal municipal waste.

#### **NOTES:**

In the event of a device-related „serious incident“ that directly or indirectly led to, could have led to, or is likely to lead to any of the following:


- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above „serious incident“ must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

#### **NOTES:**

In the event of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, a healthcare professional should be consulted.

### **DIAGRAM OF HOW TO PUT ON THE ORTHOSIS:**

### **REMOVAL**

Remove the orthosis by following all the above steps in reverse order.

### **SIZING TABLE** (To fit the correct size, measure the circumference above the ankle):

- S (19-22cm),
- M (22-26cm),
- L (26-30cm)

### **TARGET PATIENT GROUP**

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/sizes/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

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### WARRANTY CARD

Model: .....  
Serial number: .....

Point of sale stamp  
and a legible signature of the Seller: .....

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. **THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.**
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.  
03-068 Warsaw, ul. Zawiślańska 43  
Tel. 22 518 36 00, Fax 22 518 36 30  
[www.antar.net](http://www.antar.net)

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