AT53039

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USER MANUAL KNEE ORTHOSIS

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NÁVOD NA POUŽÍVANIE ORTÉZY KOLENA

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ENG

Knee brace with spiral side reinforcements and silicone patella reinforcement made of breathable fabric strengthens the joint, provides compression, stabilization and gently warms the knee joint.

CAUTION:

When a product-related "serious incident" occurs, which directly or indirectly led, could or may have led to any of the following events:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person, or
- c) a serious threat to public health

this "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION:

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of a medical device, please consult a healthcare professional.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose

ATTENTION: The manufacturer is not responsible for damage caused by neglect of maintenance, inadequate servicing or resulting from failure to comply with the recommendations contained in this manual.

INTENDED USE / INDICATIONS:

- diseases of the patellar ligament (jumper's knee)
- rehabilitation after various diseases of the knee joint
- a precautionary measure in the event of a weakened knee
- slight sprains of the knee joint

CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

FUNCTIONS:

- the orthosis ensures the functioning of the joint in a controlled, optimal, anatomical range of motion
- spiral reinforcements provide additional stabilization, joint relief, gentle massage, ensure better adhesion of the orthosis and facilitate insertion of the orthosis (no possibility of removing the reinforcements)
- silicone reinforcement ensures stabilization of the kneecap and prevents its pathological deformation

USE:

The brace is put on by pulling it over the leg.

Length of the orthosis: 30 cm

SIZES:

M: 31.7 - 36.7cm

L: 36.8 - 41.8 cm

XL: 41.9 - 47.1 cm

XXL 47.2 - 54.5 cm

CLEANING INSTRUCTIONS:

The brace should be hand washed in water with mild soap at a maximum temperature of 30 ° C. Do not use bleach, iron or dry clean. Dry in the shade, avoid direct sunlight.

WARNINGS: do not use the product if you have used warming lotions or ointments, do not use in the event of skin lesions.

In order to choose the right size, measure the circumference of the knee joint in the upright position in the middle part.

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the medical device can be disposed of as normal household waste, except for electrical products - proceed in the manner appropriate for the disposal of electrical and electronic equipment.

SHELF LIFE

Based on the tests performed, the manufacturer indicates a 5-year shelf life of the medical device.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślńska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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