AT53062/AT53063
INSTRUKCJA OBSŁUGI ORTEZA KOLANA
KNEE ORTHOSIS USER MANUAL
GEBRAUCHSANWEISUNG FÜR DIE KNIEORTHESE
NÁVOD K POUŽITÍ KOLENNÍ ORTÉZY
NÁVOD NA POUŽITIE ORTÉZY NA KOLENO
ORTHÈSE DU GENOU MODE D'EMPLOI
GEbruikSAANWIJZING KNIEORTHESE
ORTESIS DE RODILLA INSTRUCCIONES DE USO
ISTRUZIONI PER L'USO DELL'ORTESI DI GINOCCHIO
KNÄORTOS BRUKSANVISNING

PL

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EN

Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions and recommendations, please contact your health care professional or dealer to avoid product damage or personal injury.

NOTE: Inspect all parts for damage caused in transit. If you notice such damage, DO NOT use the product. Contact your dealer for more information.

ATTENTION: it is forbidden to use the product other than for its intended purpose.

NOTE: Before use, check that no part is damaged and looks correct.

NOTE: The manufacturer is not responsible for damage caused by negligent maintenance, inadequate servicing or resulting from failure to follow the recommendations in this manual.

NOTES:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person or
- (b) temporary or permanent deterioration of the health of the patient, user or other person or
- (c) a serious threat to public health

the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions, or other disturbing symptoms that are unclear to you related to the use of the medical device, consult a healthcare professional.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this product is intended (see the section on the purpose of the product of this manual). The product may be purchased by the user independently or on the recommendation of a doctor, therapist or other specialist. Whether purchasing the product on your own or as indicated by a doctor/therapist/other specialist, you should consider the available sizes/necessary functions/sizes/variants of the device, indications and contraindications for use, as well as the information provided by the manufacturer.

COMPONENTS

Orthosis, instruction manual

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling.

APPLICATION/INDICATIONS AT53062

- acute/chronic pain syndromes
- Treatment of knee joint injuries requiring limitation of range of motion
- Injuries and dysfunctions without immobilizing the joint
- support of the ligamentous apparatus of the knee joint

USES/INDICATIONS AT53063

- treatment of knee joint injuries
- Injuries and dysfunctions without immobilizing the joint
- support of the ligamentous apparatus of the knee joint

DISPOSAL OF THE PRODUCT

Dispose of the product in accordance with the applicable disposal conditions.

CLEANING/MAINTENANCE

Hand wash the product at up to 40°C with a mild detergent and then let it dry. Do not use harsh objects, solvents or strong detergents. Do not tumble dry, dry clean.
do not use bleach.


STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

SIZES AT53063 (each size has an adjustment range) it is recommended to choose a size by trying it on.

To choose the right size, measure the circumference of the thigh above the knee and the circumference of the calf below the knee.

SIZE	S	M	L	XL
Thigh circumference above the knee (cm)	39-47	47-53	53-60	60-67
Calf circumference below the knee (cm)	29-37	37-43	43-50	50-57

Height of the orthosis AT53063 - 30 cm

INSTRUCTIONS FOR USE AT53063

NOTE: The product should be worn and adjusted in the presence of a doctor or specialist on a slightly bent or straight leg.

1. Unfasten all Velcro, lay the orthosis flat with the wider side facing up.
2. Place the orthosis at knee height - so that the circular opening is at the back of the knee.
3. Then fasten all the Velcro connection to adjust the pressure in the correct order:
4. Strap 1: Fasten the orthosis on the lower thigh.
5. Strap 2: Fasten the orthosis on the upper calf.

NOTE: Fastening too tightly may lead to slowed blood circulation or a downward displacement of the bracket. Fastening too loosely will cause the orthosis to fall off.

FEATURES AND FUNCTIONS AT53063


- Made of high quality materials, which ensures high comfort of use
- Supports stabilization of the knee joint and patella

It has a wide range of adjustment thanks to the Velcro fastening system

SIZE AT53062 (each size has a range of adjustment) it is recommended to choose a size by trying it on.

To choose the right size, measure the circumference of the thigh above the knee and the circumference of the calf below the knee.

SIZE	S	M	L
Thigh circumference above the knee (cm)	39-47	47-53	53-60
Calf circumference below the knee (cm)	29-37	37-43	43-50

MOBILITY ADJUSTMENT RANGE AT53062

- Bending settings: 0, 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 110, 120°

- Upright settings: 0, 10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 110, 120°

Height of the orthosis AT53062 - 34 cm

INSTRUCTIONS FOR USE AT53062

NOTE: The product should be worn and adjusted in the presence of a doctor or specialist on a slightly bent or straight leg.

1. Unfasten all Velcro, lay the orthosis flat.
 2. Place the orthosis at knee height - so that the clocks are at knee joint height on the left and right sides.
 3. Fasten the orthosis in front with a Velcro connection.
 4. Then fasten all 4 Velcro straps to adjust the pressure in the correct order:
Strap 1: Fasten the strap on the lower thigh.
Strap 2: Fasten the strap on the upper calf.
Strap 3: Fasten the strap on the upper thigh.
Strap 4: Fasten the strap on the lower part of the calf. Fastening too tightly can lead to slowed blood circulation or shifting of the bracket downward. Fastening too loosely will cause the orthotic to fall off.
- Note:** Viewed from the side, the hinge should be slightly behind the midline of the knee.
5. The orthosis can be used in a locked position of the knee or with the ability to bend. To adjust the range of flexion angle adjustment, unlock the adjustability with the button (fig.1)
 6. Then set the desired values on the timer (starting from 0) and lock the timer by moving the button. (fig.2)

1.

2.

IMPORTANT:

- 1) The product is intended for use by one person.
- 2) Both hinges must be set to the same upright and/or flexion settings.
- 3) The patient should consult a physician before adjusting the orthosis.
- 4) Proper rehabilitation and activity improvement are also essential processes in safe treatment planning.

FEATURES AND FUNCTIONS AT53062

- Ability to stabilize and lock the knee joint at a specific angle
- Orthosis has adjustable range of knee joint flexion angle

It has a wide range of adjustment thanks to the Velcro fastening system

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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