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EN

Thank you for purchasing our product. Please carefully read the user manual. Do not use the product before reviewing and understanding these instructions. If you do not understand the warnings, notices, or recommendations, contact a healthcare professional or the seller to avoid product damage or personal injury.

WARNING:

Check all components for damage caused during transport. If such damage is observed, DO NOT use the product. For more information, contact the seller.

TARGET PATIENT GROUP

This product is intended for individuals dealing with diseases, dysfunctions, or injuries for which this product is designed to treat/rehabilitate or compensate (refer to the purpose of the product section in this manual). The product may be purchased by the user independently or on the recommendation of a doctor, therapist, or other specialist. Regardless of whether the product is purchased independently or recommended by a healthcare professional, the available sizes, required features, product variants, indications and contraindications for use, as well as information provided by the manufacturer, must be considered.

APPLICATION/INDICATIONS

The orthosis is a medical device designed to stabilize, unload, or correct the spine, used to treat injuries, diseases, or dysfunctions of the musculoskeletal system within the spine.

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairments) preventing safe use of the product.

Skin damage at the site of orthosis application: open wounds, abrasions, skin infections, pressure sores.

Severe cardiovascular diseases: impeded blood flow caused by orthosis pressure. Diseases requiring avoidance of chest and abdominal compression (e.g., abdominal aortic aneurysm).

Improper orthosis fitting: An orthosis that is not appropriately adjusted or correctly worn may worsen dysfunctions instead of alleviating them.

Neurological contraindications: Neurological conditions where stabilization may exacerbate symptoms (e.g., spasticity or neuropathic pain).

Osteoporosis or advanced bone fragility: In some cases, the orthosis may impose undesirable stress on fragile bone structures.

Respiratory diseases: Difficulty breathing, especially when using orthoses covering the chest (e.g., asthma, chronic obstructive pulmonary disease).

Inflammation or swelling at the orthosis application site: Chronic inflammation that may worsen due to orthosis use.

Lack of patient consent or hygiene maintenance difficulties: Improper skin care under the orthosis may lead to complications.

Individual anatomical restrictions: Spinal or chest deformities preventing proper orthosis fitting.

FEATURES AND SIZE

Material:

Main material: 80% Polyester, 20% Nylon

Metal support: Steel and aluminum alloy plates

Velcro: 100% Nylon

Straps: 70% Polyester, 30% Latex

Easy to use

STABILIZING ELEMENTS:

There are three non-adjustable plates that can be removed from the pockets.

Material of the large plate in the middle: Aluminum alloy

Material of the small plate: Steel plate.

SIZES

Length: M - 90 cm / L - 100 cm / XL - 110 cm / 2XL - 120 cm / 3XL - 130 cm

Waist: M - 70-85 cm / L - 85-95 cm / XL - 90-105 cm / 2XL - 100-115 cm / 3XL - 110-125 cm

COMPONENTS

Spinal orthosis, User manual

WARNING:

In the event of a "serious incident" related to the product that directly or indirectly caused, could cause, or may cause any of the following:

a) death of a patient, user, or another person, or
 b) temporary or permanent deterioration of health of the patient, user, or another person, or
 c) serious threat to public health,
 the above "serious incident" must be reported to the manufacturer and the competent authority of the member state where the user or patient resides. In Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

WARNING:

If pain, allergic reactions, or other disturbing symptoms occur during the use of the medical device, consult a healthcare professional.

WARNING:

The product must not be used in a manner inconsistent with its intended purpose.

USAGE

Place the orthosis on the spine, adjust the compression using the straps, and fasten it with Velcro.

Rear Pelotte

Please familiarize yourself with the position of the rear pelotte in the attached image. The rear pelotte should adhere to the material in the center of the strap.

Material: EVA and polyester fabric.

CLEANING

The product should be hand-washed with mild detergent at temperatures up to 40°C. Do not machine wash. Avoid using strong detergents, bleach, or abrasive materials. Dry flat in the shade; do not expose to direct sunlight. Do not tumble dry. Do not dry-clean. Do not iron.

STORAGE

Store the product in a dry place, away from direct sunlight.

DISPOSAL

When no longer in use, the product should be disposed of as regular municipal waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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