AT53108/AT53109

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EN

Thank you for purchasing our product. Please read this user manual carefully before use. The medical device must not be used before you have fully read and understood these instructions. If you do not understand any of the warnings, notes, or recommendations, please contact a healthcare professional or the distributor to avoid product damage or personal injury. NOTICE Check all parts for possible transport damage before use. If you notice any damage, DO NOT use the product. For more information, contact your distributor.

INTENDED PATIENT GROUP: This medical device is intended for individuals suffering from diseases, dysfunctions, or injuries for which the device is designed (see the "Intended Use" section of this manual). The product may be purchased independently by the user or based on the recommendation of a physician, therapist, or other healthcare professional. In both cases — whether purchased directly or recommended by a specialist — it is important to consider the available sizes, required functions, product variants, as well as the indications, contraindications, and information provided by the manufacturer.

INTENDED USE: The pneumatic ankle orthosis is designed to immobilize, stabilize, and support the lower limb — particularly the calf, ankle, and foot — in various pathological and post-traumatic conditions.

Applications include: Post-injury stabilization: Recommended for fractures, sprains, or other injuries of the foot and/or lower leg where immobilization and support are required during recovery. Post-operative recovery: Used after surgical procedures involving the lower limb to ensure stability and protection of the operated area. Chronic conditions: Suitable for treating chronic ankle instability, foot deformities, or other conditions requiring structural support.

Indications: Stable foot or ankle fractures; severe ankle sprains; post-operative lower limb recovery.

Contraindications: Physical or mental limitations (e.g., visual impairment) preventing safe use of the product. Paresthesia (sensory disturbances). Advanced varicose veins. Skin diseases or injuries in the area of application. Lymphatic drainage disorders.

PACKAGE CONTENTS: Pneumatic ankle orthosis, User manual

FEATURES

Dual air cushions on both sides independently control the compression level as needed.

Air pump allows efficient inflation and deflation.

Anatomically shaped outer boot design provides excellent stability.

Special contoured sole allows natural, comfortable walking.

SIZES (shoe size in cm): XS: <21.5 S: 22–25.5 M: 26–28 L: 28–30.5 XL: >31 Height of AT53108: 43 cm Height of AT53109: 28 cm

MATERIAL COMPOSITION

Outer boot structure: PP 56%, EVA 13%, PA 10%, POM 3%, steel 2%, TPU 10%, PVC 6% Inner liner: PA 30%, PU 60%, polyester 10%

WARNING: Do not wear the orthosis continuously for more than 24 hours without a break.

INSTRUCTIONS FOR USE

Insert the foot and lower leg into the foam liner and fasten securely. If necessary, use additional pads to fill gaps. Place the wrapped foot into the orthosis and attach the front panel for better protection and support. Fasten the straps at the ankle, foot, and calf in the order shown in the illustration below. Tighten the straps until comfortably snug.

AIR PRESSURE ADJUSTMENT

Press the air pump to inflate. Press the valve to release air.

WARNING

In the event of a serious incident related to the medical device that has led, might have led, or could lead to any of the following:

- a) the death of a patient, user, or other person;
- b) a temporary or permanent serious deterioration of the health of a patient, user, or other person;
- c) a serious threat to public health;

the serious incident must be reported to the manufacturer and to the competent authority of the EU Member State in which the user or patient resides. CAUTION: If pain, allergic reactions, or other concerning or unclear symptoms occur during use of the medical device, consult a healthcare professional immediately. It is forbidden to use the product in any way other than its intended purpose.

CLEANING

Fasten all straps before washing. Hand wash in cold water with mild soap. Do not dry clean, bleach, or tumble dry. Allow to air dry naturally — avoid heat sources.

STORAGE

Store the product in a dry place, away from direct sunlight. Temperature: -5°C to 40°C Humidity: 35%–75%

DISPOSAL

After end of use, dispose of the product as general household waste in compliance with local waste segregation regulations.

Be environmentally responsible — recycle when possible.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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