
700214-700215-700216

Mode d'emploi
Bedienungsanleitung
User Manual
Instrucciones de uso
Upute za uporabu
Instruções de utilização
Kullanım kılavuzu

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CONTENTS

1. PURPOSE AND REUSE OF THE CHAIR

1.1 Purpose, indications and counter-indications

2. SECURITY INFORMATION

2.1 Pre-use checks

2.2 Indoor use

2.3 Outdoor use

2.4 Transporting the chair

2.5 Risk of pinching

2.6 Spare parts

3. RECEPTION AND DESCRIPTION

3.1 Unpacking and contents of the delivery

3.2 Product description

3.3 Folding and unfolding

3.4 Installing the footrests

3.4.1 Fitting the footrest

3.4.2 Removing the footrest

3.4.3 Height adjustment of the footrest

3.5 Fitting the armrests

3.5.1 Replacing the armrests

3.6 using the brake

3.6.1 Parking brake

3.6.2 Footbrake

3.6.3 Brake adjustment

3.7 Wheels

4. PRACTICAL ADVICE FOR A GOOD USE

4.1 Moving and braking

4.2 Getting in and out of the chair

4.2.1 Getting in and out of the chair from the front

4.2.2 Lateral transfer

4.3 Overcoming obstacles, ascending or descending slopes

4.4 Lifting the chair, climbing stairs

4.5 Picking up an object

4.6 Stability when moving or stationary

5. DAILY CARE AND ADJUSTMENT INSTRUCTIONS

5.1 Maintenance instructions before each use

5.1.1 Braking system

5.1.2 Tyres

5.1.3 General condition

5.2 Maintenance, adjustment and regular checks

5.3 Conditions of storage and use

5.4 Cleaning

6. DISPOSAL

6.1 Checks and repairs

6.2 Storage

6.3 End of life of the chair

7. REUSE

8. GUARANTEE

9. TECHNICAL DATA

9.1 Essential characteristics of the chairs

9.2 Dimensions

9.3 Recommendations

9.4. Maintenance booklet

10. GUARANTEE VOUCHER

CONTENTS

By reading this guide beforehand, you will be familiar with the features and recommendations that are essential for the optimal use and maintenance of your device. Do not use this device until you have read this manual for the first time. Always carry out a general inspection of your device before use. If you are unsure of any aspect of your device, contact your nearest dealer. As an authorised professional, he will be able to advise you and check your chair before and after use.

GENERAL SAFETY INSTRUCTIONS

- Ask your specialist dealer or medical personnel to teach you how to use and maintain this device.
- Do not store outside and do not leave in the sun (Caution: do not leave in the car at high temperatures!).
- Do not use at extreme temperatures ($> 38^{\circ}\text{C}$ / 100°F or $< 0^{\circ}\text{C}$ / 32°F).
- Regularly check all mechanical parts.
- Regularly check all screw connections for tightness.
- VDo not use for other purposes than transfer chair use: e.g. do not stand on the transfer chair!
- Observe the maximum weight.
- No changes may be made to the product. Only use original spare parts and have repairs carried out by qualified personnel.
- Once you have used the transfer chair, make sure that it is locked and stable in the open position.
- Do not attach anything to your transfer chair.
- If you have not used the transfer chair for a long time, have it checked by qualified personnel first.
- Check that the transfer chair is in perfect condition before each use.
- If possible, only use the transfer chair on a horizontal, flat and stable floor. Transfer chairs have only limited resistance to tipping and sliding.
- Sitting on the armrests is strictly forbidden.
- Do not bend over to avoid falling.
- This product may only be adjusted by qualified personnel. In case of doubt, contact your specialist dealer.
- **In the event of prolonged use, a pressure sore cushion must be used on the chair to prevent the development of pressure sores.**
- The manufacturer and the competent authority of the Member State in which the user and/or patient is established must be notified of any serious incidents that occur in connection with the device.

1. PURPOSE OF USE, INDICATIONS AND CONTRAINDICATIONS

1.1 Purpose of use

Our Aliv'up Plus transfer chair is meant to accompany and move from one place to another people who don't have specific postural problems and whose weight don't exceed 115kg. This chair is not conceived for very young children or patients whose weight is inferior to 40kg or for other uses than the transport of a person. This chair needs the presence of a third person, responsible and physically able to move the chair.

The Aliv'up Plus chair is a light chair, manoeuvrable and usable on a plain surface indoors and outdoors. Foldable, its footprint is reduced to be able to be transported and be stored easily.

Only healthcare professionals are authorised to prescribe and recommend transfer chairs for use in accordance with the patient's needs.

Adaptation and adjustment of the transfer chair must be carried out by an authorised professional.

This approval is given by the Cerah in Woippy 57 for France.

Indications for use :

Inability to walk or severely limited ability to walk as part of the basic need to move around in one's own home and to leave home for a short walk in the open air or to reach places usually close to home where routine business has to be carried out.

ALIV'UP PLUS is a lightweight transfer chair and is suitable when disabled people are unable to move by themselves, without assistance, for example to climb stairs, to enter in a plane or a coach.

Contra-indications :

Weakness of the trunk	Postural deficiency	Herniated disc	Lower back pain
Ischial pressure sores	Fractures of vertebrae	Visual deficiency	Uncontrolled movements
Cognitive impairment	Some pelvic fractures	Surgery or post-operative pelvic pain	Problems with the femurs

2. GENERAL USE

2.1 Pre-use checks

Before each use, the user or accompanying person must check that the wheels, tyres and bearings are in good condition and that the brakes allow the transfer chair to be stopped effectively.

Never use the transfer chair if a clamping element is unscrewed, if a wheel or part of the transfer chair has abnormal play or if the upholstery is torn or loose. Never use the chair if it is not in good condition.

2.2 Indoor use

When moving the transfer chair, it is important to take into account its own dimensions, as well as any objects, steps and obstacles, so as to avoid hindering or injuring the patient's movements.

2.3 Outdoor use

When moving the transfer chair, it is important to take into account its own dimensions, as well as any objects, steps and obstacles, so as to avoid hindering or injuring the patient's movements.

2.4 Transporting the transfer chair

It is not advisable to lift the chair with the patient sitting on it. To transport the chair, fold the frame and lift it using stable structural elements (armrests, backrest handles, frame) excluding footrests or wheels.

When transporting the chair in a vehicle, it is advisable to remove the footrests to reduce the space required for storage in a boot. In no way should the transfer chair be used as a substitute for the chairs in the vehicle, the patient must be transferred from the transfer chair to the vehicle by first locking the transfer chair brakes.

2.5 Risk of pinching

When installing or storing the transfer chair, care must be taken to avoid possible pinching due to the adjustment of the footrests, the removal of the wheels, the folding of the seat....Ask for advice from your authorised dealer, who will explain how to handle your transfer chair safely.

2.6 Spare parts

The replacement of defective or altered components should only be done with parts guaranteed by the manufacturer to a professional. Any replacement with non-compliant parts may entail a risk for the patient and nullify the manufacturer's warranty.

3. RECEIPT & DESCRIPTION

3.1 Unpacking and contents of the delivery

The chairs are checked at various stages of their manufacture and then packed in a special cardboard packaging. On receipt, check that the cardboard box has not been damaged in any way, then when unpacking, check that the accessories are present, such as the instructions for use and the guarantee booklet. These documents must be given to the user so that he or she is aware of the conditions of use and warranty of the chair.

If in doubt about its operation, contact your dealer for advice, and repair if necessary. A summary table at the end of the manual contains essential information about the chairs and their limits of use.

3.2 Product description

Our Aliv'Up Plus are controlled at different steps of their production then packed in a cardboard packaging. Upon receipt, please check that the packaging is not damaged. When unpacking the product, make sure that the instruction manual and the warranty booklet are present. These documents are to be delivered to the user imperatively for him or her to be aware of the use conditions and the warranty of the product. The chair is equipped with :- A seat with a backrest foldable half-height, a canva to sit on it, two armrests fixed laterally and a fastening strap.

- A foldable chassis with two lateral flanks and a cross system that allows to fold and unfold the chair.
- Two multidirectionnal front wheels
- Two back wheels that can be blocked with the brakes located on the sides and/or at the top of the backrest.
- Two adjustable in height and removable footrests fixed on each flank of the chair.

After having checked the conformity of the chair, the certified professional will adjust the footrests according to the user's height, as well as the strap that is used as a support belt. It is the professional's responsibility to inform and to make sure that the user understood well the safety or use instructions and warnings described in this instruction manual. Environment : Recyclable, the cardboard packaging must be taken to a sorting center.

It is the responsibility of the professional to inform and ensure that the user understands the safety or operating instructions and warnings described in this manual.


3.3 Folding and unfolding

Beware of the risk of pinching when folding and unfolding the chair.

Folding and Unfolding :

To fold the transfer chair : stand behind the transfer chair, then lift up the two handles located on the two sides of the seat (1) and (2).

To unfold the transfer chair: stand behind the transfer chair, then move the two sides of the transfer chair (2) apart and finally press down on each side of the canvas seat (3) so that it is in contact with the 3 supports on each side of the transfer chair. Footrests can be folded on side or not during this operation. Straighten the backrest, locking the upper part with the two brackets (see next step «Folding the backrest»).

Folding the backrest:

To fold the backrest, press on the lever (2) located on both sides of the frame and fold the backrest (3).


3.4 Installation of footrests

Both footrests are delivered with the paddles mounted. Both footrests can be folded up in the sitting position.

3.4.1 Fitting the footrest

1. Hold the footrests on the sides of the transfer chair frame.
2. Insert parts [A] and [B] of the footrest into brackets [C] and [D] of the frame.
3. Swivel the footrest inwards until it audibly locks into its final position.
4. Return the footrest paddles to the horizontal position.


5. Replace the footbelt


3.4.2 Removal of the footrest

1. Push lever [E] to unlock the footrest.
2. Swivel the footrest to the outside of the transfer chair by 90°.
3. Pull the footrest upwards out of the tube to remove the footrest from brackets [C] and [D].


3.4.3 Footrest height adjustment, see [A]

Caution : Risk of accident if the footrests are set too low.

1. Unscrew the screws [F] at the footrests with the hexagon key provided.
2. Pull the tube downwards and adjust it to the desired height.
3. Tighten the screws again. Make sure that both paddles have the same height adjustment.

3.5 Fitting the armrests

The armrests can be tilted backwards.

3.5.1 Replacing the armrests

The armrests are easy to replace, simply unscrew and remove the nut [I], then adjust the new armrest and screw it back on again

3.6 Using the brake

The brakes can be operated by means of the wheel lock lever or the parking brake.


Caution: The brakes are only designed to stop the transfer chair and not to brake the transfer chair in motion. Checks and adjustments must be carried out regularly by an authorised dealer to ensure maximum safety for the user. When the transfer chair is on a sloping surface, release the brakes one by one.

3.6.1 Parking brakes

CAUTION: Make sure that the transfer chair is on a flat, horizontal surface before releasing the brakes. Never release both brakes simultaneously.

Parking brakes (1) and (2):

1. To break, push the lever downwards
2. Release the brake by pulling lever up

3.6.2 Foot brake

1. To break, push on the lever downwards (1).
2. Release the brake by pulling the lever upwards (2).

3.6.3 Brake adjustment

The brake pressure on the tyre can be adjusted by unscrewing the fixing screws and sliding the brake on the tube..

3.7 Wheels

The rear wheels (300 mm) and the front wheels (200mm) are delivered mounted. The wheels are fitted with puncture-proof PU tyres. Check that it is not damaged before each use.

4. PRACTICAL ADVICE FOR A GOOD USE

4.1 Moving and entering/exiting the chair

The movement is carried out with the help of an assistant. When driving the transfer chair, assistant must be careful.

4.2 Getting in and out of the transfer chair

When sitting in or getting out of the chair, make sure that it is placed on a flat, horizontal surface without obstacles to avoid any risk of falling.

To get into the chair, make sure the patient is placed close to the chair, block the rear wheels, the parking brakes and raise the footrest paddles to avoid hitting them during transfer. It is also possible to lift the retractable armrests for easier access to the seat or for a side transfer.

To get out of the chair, lift the pallets and the armrests with the wheels locked.


If necessary, call for outside help to assist you during the transfer.

4.2.1 Entering and exiting from the front

1. Fold the footrests upwards and swivel the leg rests to the side.
2. bring the chair and the user as close as possible to each other from the front.
3. Brake the rear wheels and parking brake.
4. Slide on the seat and sit backwards.
5. Return the footrest paddles to the horizontal position.

4.2.2 Lateral transfer

1. Fold up the footrests and armrests and remove the legrests.
2. Bring the chair as close as possible to the side of the user
3. Brake the rear wheels.
4. Slide on the seat and sit backwards.
5. Return the footrest paddles to the horizontal position.


Caution: It is strictly forbidden to lean on the footrest paddles when trying to get out of the chair. Danger of tipping over and injury!

4.3 Overcoming obstacles, uphill or downhill slopes

When driving on a slope, please note that there is already a risk of tipping over on slopes with a gradient of 8% (5°), and leaning out of the chair increases the risk of tipping over. The lateral inclination should not exceed 15°.

For safety reasons, the assistant should drive the chair over obstacles such as pavements and steps very carefully or with the help of a second accompanying person. The accompanying person must be able to brake the chair on his or her own.

Depending on the effort or type of obstacle, the manoeuvre must be performed forwards or backwards.

4.4 Lifting the chair, climbing stairs

Never lift the chair using the armrests, but use the horizontal tubes on the sides of the seat. The help of two caregivers is essential for the user's safety.

When climbing stairs, it is advisable, as when descending a slope, to reverse the position of the transfer chair so that each caregiver is positioned on either side of the transfer chair. The caregiver at the lowest point will grab the lower front part of the transfer chair (never the footrests), the caregiver at the top of the stairs will use the handles to tip the transfer chair towards him or her to allow the user to climb the stairs.

4.5 Pick up an object

When searching for an object in front of the chair, the multidirectional castors should be positioned forwards to increase the stability of the chair, for an object behind, the castors should be positioned backwards. Under no circumstances should the chair be tipped on the rear wheels without the assistance of a third party. In any case, it is recommended to lock the brakes on the rear wheels beforehand.

4.6 Stability in motion or at rest

Designed to remain stable when moving around and when stationary, the stability of the chair is nevertheless dependent on maintaining the patient's centre of gravity within an area between the four points where the wheels rest on the floor. Both the patient and the caregiver must ensure that this is maintained even when reaching an object. Only the arms may be extended beyond the chair. The buttocks must always rest on the seat and under no circumstances should they be supported on the paddles.

When the chair is unfolded, the seat and backrest canvases should be stretched. Under no circumstances should the chair be used if they are cut, torn or in danger of breaking. If in doubt, always replace them.

EN 5. DAILY CARE AND ADJUSTMENT INSTRUCTIONS

5.1 Maintenance before each use

5.1.1 Braking system

As an essential safety feature, the brakes must be checked before each use to validate their effectiveness in keeping the transfer chair motionless. Adjustment is carried out by loosening/tightening the two screws allowing the flange to slide along the tube. The braking system must maintain the expected efficiency. In case of doubt, any parts found to be faulty should be replaced.

5.1.2. Tyres

Tyres must be in good condition, free of tears or tearing, and must turn easily without noise or squealing.

5.1.3 General condition

Wheels : The wheels and their axles must be checked for tyre quality, deformation and axle play. The wheels must turn easily, without any particular noise and perpendicularly to the axis on flat ground, and stop immediately using the brake or in contact with the ground (without pushing the transfer chair). The rear and front wheels do not need to be inflated on chairs with full wheels.

Footrests : the footrest paddles must not have any cracks or be broken, and must be able to articulate correctly on their axis of rotation. Once installed, the footrest locking clips on each side must keep the footrests firmly attached to the frame. The footrest paddles must be changed in the event of cracks or malfunction. The footrest height adjustment tube must be able to slide unhindered inside the top tube and lock securely. The height of the footrests can be adjusted using the clamping screw.

Screws and bolts : All screws and bolts as well as the clamping parts must be checked and replaced if necessary. The threads and screw cross-sections must be correct and free of wear and rust.

5.2. Maintenance table

For the perfect safety and longevity of the chair, it is recommended to have the chair serviced by an approved professional at least once a year, or as soon as a problem is encountered.

Please note that only original transfer chair components should be used to repair the transfer chair, otherwise the manufacturer's warranty will be invalidated.

A complete test of the transfer chair must be carried out before returning it to the user, and the intervention must be recorded in the maintenance and warranty booklet.

Operations	Frequency	User	Authorised distributor (Cerah)
Installation and commissioning	1st use		commissioning and training of the user and the attendant
Component replacement			to be carried out by the authorised distributor
Wheel tyre check	before each use	mandatory	in case of doubt, contact your authorised distributor
Brake efficiency	before each use	mandatory	in case of doubt, contact your authorised distributor
Visual inspection, tightening and cleaning	once a month or when necessary	highly recommended	in case of doubt, contact your authorised distributor
Revision	annual		to be carried out by the authorised distributor
Repair	as required		to be carried out by the authorised distributor
Supply of spare parts	as required		to be carried out by the authorised distributor

5.3 Conditions of storage and use

	Temperature		Moisture	
	mini	maxi	mini	maxi
Storage	5° C	+ 50° C	0%	80%
Use	0° C	+ 40° C	0%	80%

5.4 Cleaning

For a safe use and longer longevity, the Aliv'Up chair must be regularly cleaned up. Remove all traces of dirt on the chassis, the brakes and the rotation axis with the help of a soft cloth, and non corrosive household product (liquid household cleaner, soapy water) before rinsing it. The seat, backrest and armrests as well as the handles situated at the back must be cleaned up too for a better comfort of use.

When changing users, it is recommended to disinfect the chair more thoroughly with disinfectant solutions such as Surfanios[®], following the product's instructions.

6. DISPOSAL

Do not throw away the chair or its accessories when they are no longer usable, return them to an approved waste disposal centre. Ask your specialist dealer about regional disposal regulations.

7. REUSE

If used by different people, ALIV'UP PLUS chairs must be inspected for servicing, maintenance and possible adjustment for the new user by an authorised professional. The user manual and the original documents and labels must accompany the chair and be handed over to the new user.

8. WARRANTY CONDITIONS

Transfer chairs are insured for 2 years against any manufacturing defect under normal use. In the event of repair or replacement of parts, the guarantee will be maintained subject to the intervention of an approved professional and subject to the use of components supplied by Herdegen.

The warranty does not cover normal wear and tear of materials and components.

The manufacturer's guarantee excludes any damage and any direct or indirect implication.

The transfer chairs have been tested internally and meet the requirements of MDR 2017/745.

The warranty booklet is supplied with each transfer chair and must be kept and completed at each intervention.

The transfer chairs have been validated by internal tests and meet the current specifications established by the Ministry of Health and applied by the CERAH certified body.

9. CHARACTERISTICS AND RECOMMENDATIONS

9.1 Main characteristics of the transfer chairs

Brakes : manual

Surface finish : epoxy and chrome painting for footrest adjustment tubes


Seat, backrest and armrests : Non-fire vinyle

Aluminium frame.

Polyurethane wheels without inner tube.

9.2 Dimensions

Codes	Items	A	B	C	D	E	F	G	ø wheels	Weight
700214	ALIV UP PLUS 39,5 cm	46 cm	48.5 cm	48 cm	39 x 39 cm	108 cm	75 x 27 x 64 cm	96 cm	15/30.2 cm	11.80 kg
700215	ALIV UP PLUS 43 cm		52 cm	51.5 cm	43 x 39 cm					11.94 kg
700216	ALIV UP PLUS 48 cm		57 cm	56.5 cm	48 x 39 cm					12.04 kg

A= Height seat-to-floor
B= Overall width
C= Width between armrests
D= Seat dimensions
E= Depth
F= Dimensions folded
G= Height handles-to-floor

9.3 Recommendations

User weight: min. 40 kg, max. 115 kg. Duration of use: 5 years.

A doctor or an authorised medical professional can recommend and prescribe the model according to the pathology and needs of the patient, taking into account the number mentioned above.

EN 9.4. Maintenance booklet

date and stamp , description	date and stamp , description	date and stamp , description
date and stamp , description	date and stamp , description	date and stamp , description
date and stamp , description	date and stamp , description	date and stamp , description

10. Warranty Voucher

Warranty Voucher	
Product :	Eagle
Reference:	_____
Serial number *:	_____
* (to be fill in by the seller)	
Retailer:	
Date & Signature:	

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HERDEGEN SAS