
PNEUMATIK ANTI DECUBITUS AIR CELL CUSHION

Models: AT52108, AT52109, AT52110

Pneumatic anti-decubitus cushion

The pneumatic anti-decubitus cushion allows the user a comfortable and stable sitting position. The Single Room Pneumatic Pressure Relief Pillow is designed to evenly distribute the patient's weight and adjust the pillow to the user's current body position. The pillow is soft, lightweight, easy to carry and clean.

The product allows the sacrum to be fully immersed in the pillow, reducing pressure on the coccyx and sciatica. The cushion prevents pressure points, prevents pain during movement, reduces the risk of damage to skin tissue and optimizes tissue relief.

- Combined air chambers are easy to inflate and adjust.
- A specially designed structure of air-filled spaces allows air to flow between them, optimizing the distribution of body pressure on the pad.
- Each pillow contains a valve that is used to adjust the amount of air in the pillow and therefore the softness of the pillow.

Indications:

- Patients with hip ulcers.
- Patients with pelvic rotation or pelvic changes.
- Patients with stents.
- People in wheelchairs.

Maximum user weight AT52109: 120 kg. Maximum user weight AT52108, AT52110: 136 kg

WARNING:

1. Read the instruction manual carefully before using the pillow.
2. Be careful of what the bottom part of the pillow contacts as there may be some shifting phenomenon. Be careful.
3. Do not use sharp tools when touching the cushion.
4. Inflate the cushion with the correct amount of air before use.
5. Keep the pillow away from fire and moisture.
6. Do not bend the pillow when filling air.
7. If any part of the body is left on the cushion for too long, air will escape from that area, reducing the therapeutic effect.
8. Check the cushion daily (pressure regulation)
9. The product should be stored in a dry environment.

CAPACITY

1. The cushion provides comfort and stability by evenly distributing the pressure from the hips so that the user sinks into it when sitting on this cushion.
2. The cushion is intended for contact use, particularly for patients at high risk of bedsores.
3. Filled gaps make your posture more balanced and stable.
4. The cushion is suitable for people with disabilities, the elderly, office workers, drivers, people who sit for a long time, etc.

PRODUCT SPECIFICATION:



Product code	Dimensions of the pillow Width Length Height	Product features
AT52108	44 x 38 x 5 cm	- small air spaces with no room for growth - 9 rows with 8 air spaces - Dimensions of the air wings: width 4 cm, length 4 cm - Net Weight: 0.33 kg - Material: PVC
AT52109	45 x 39 x 4 cm	- Airspaces at low altitude - Greater stability due to larger and lower spaces - 7 rows with 6 air spaces - Dimensions of the air wings: width 6 cm, length 6 cm - Net weight: 0.255 kg - Material: PVC
AT52110	45 x 39 x 6 cm	- relatively high airspaces - 7 rows with 6 air spaces - More pleasant feel due to larger and higher air spaces - Dimensions of the air wings: width 6 cm, length 6 cm - Net weight: 0.315 kg - Material: TPU

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

Air Wings: Made of PVC and TPU.


1. Ventilation: Each pillow contains a valve that can be used to adjust the amount of air in the pillow and thus the softness of the pillow.
2. Instruction Manual: Contains instructions for adjustment, maintenance and repair.


1. Fill the air bag orally until the cells are smooth.

2. After the cushion is inflated, the user can sit in a normal sitting position with correct arm and foot position.

3. Place a hand between the cushion and the user's tailbone to adjust the cushion into position.

4. Gently unscrew the valve by hand to deflate and adjust the pillow to the user.

WARNING:

In the event of a „serious incident“ related to the product that has resulted, could have resulted, or is likely to result, directly or indirectly, in any of the following situations:

- a) the death of the patient, user or other person, or
- b) temporary or permanent deterioration in the health of the patient, user or other person, or
- c) a serious risk to public health

this „serious incident“ must be reported to the manufacturer and to the competent authority of the Member State where the user or patient is present. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

WARNING:

In case of pain, allergic reactions or other disturbing and unclear symptoms related to the use of the drug, consult your doctor.

Indications

Hip ulcers and decubitus ulcers, rotation and other changes in the pelvis, hip arthroplasty, wheelchair patients.

Utilisation

Please dispose of municipal waste in accordance with the applicable disposal regulations.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

WARRANTY CARD

Model:

Serial number:

Point of sale stamp
and a legible signature of the Seller

1) ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.

2) During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.

3) ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.

4) If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.

5) Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.

6) The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.

7) The product delivered to the service in connection with the complaint cannot cause an epidemic.

8) THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.

9) Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Phone +48 22 518 36 00, Fax +48 22 518 36 30
www.antar.net e-mail: antar@antar.net

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