INSTRUKCJA OBSŁUGI TYLNY STEROWNIK
OPERATING MANUAL REAR CONTROL
BEDIENUNGSANLEITUNG HECKSTEUERUNG
NÁVOD K OBSLUZE - OVLADAČ PRO DOPROVOD
NÁVOD NA OBSLUHU ZADNÉ OVLÁDANIE
MANUEL D'UTILISATION COMMANDE ARRIERE
BEDIENINGSHANDLEIDING ACHTERBEDIENING
MANUAL DE FUNCIONAMIENTO MANDO TRASERO
MANUALE D'USO COMANDI POSTERIORI
BRUKSANVISNING BAKRE KONTROLL

IU_AT52304-1

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REAR COTROLLER OPERATING INSTRUCTIONS

1. INSTALLATION TIPS

- First install the fastener then the rear controller on wheelchair right handle using a Allen key.
- Insert the XLR plug/USB plug of the rear controller connecting line into charging port (XLR or USB port) of the main wheelchair controller and fasten the cable with a cable fixing.

2. SWITCH REAR CONTROLLER ON/OFF

- 1) Press the „POWER“ button for 2s until the rear controller is switched on.
- 2) Press the „MODE“ button to change the speed (5 speed modes).

3. USE

After steps 1 and 2 above, the rear controller operator holds the wheelchair handles with both hands and uses the right thumb to control the rotation lever. By moving the lever forward, the wheelchair moves forward and vice versa; when the lever naturally returns to the centre position, the electric wheelchair will not move.

4. PROBLEMS

Common problems

1. the wheelchair doesn't move whether the lever is moved forwards or backwards;
2. the wheelchair continues to move even when the lever is moved back to the centre position;

If these cases occur, verify the parameters of the rear controller:

1. disconnect the connecting line from the main controller;
2. hold the „Mode“ button, in the meantime connect the rear controller to the main controller; when the rear controller enters the verification mode, release the „Mode“ button and move the lever forward and backward at least twice (note: the lever must be moved to the furthest sides). When the lever returns to the central position, press the „Mode“ button again to save the data verification and complete the process.

CAUTION:

When a product-related „serious incident“ occurs, which directly or indirectly led, could or may have led to any of the following events:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in

the health of a patient, user or other person, or c) a serious threat to public health this „serious incident“ must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicines, Medical Devices and Biocidal Products.

CAUTION:

In the event of pain, allergic reactions or other worrying symptoms, unclear for the user, related to the use of a medical device, please consult a healthcare professional.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose

Destiny:

Additional driver designed for the assistant of the wheelchair user

Contraindications

Physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

TARGET PATIENT GROUP

The healthcare professional shall indicate the use of the device to adults and children on his/her own responsibility, taking into account the available variants/sizes/necessary functions/sizes and indications, bearing in mind the information provided by the manufacturer.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the medical device can be disposed of as normal household waste, except for electrical products - proceed in the manner appropriate for the disposal of electrical and electronic equipment.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 **ANTAR Sp. J.**
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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