AT53006, AT53030

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ENG

GENERAL WARNING

Do not use the product before reading and understanding these instructions. If you do not understand the warnings, cautions or instructions, contact your healthcare professional or dealer to avoid damage to the product.

PURPOSE

Single-chamber plastic shin and foot orthosis with pneumatic stabilisation. It is used for fractures and sprains in the foot and also in the case of damage to the Achilles tendon instead of a plaster bandage.

INDICATIONS


- Post-operative conditions of the foot or ankle
- Stabilisation of foot or ankle fractures
- Moderate to severe ankle sprains

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling.


USE

- 1 Place the foot and leg in the foam insole and tighten so that the insole is snug against the foot.
2. for a better fit, place cushions inside the insole or between the insole and the casing.
Place the foot wrapped in the insole in the shoe and pull the vertical supports on both sides slightly outwards. The lateral axis of the shin should be aligned on both sides to the supports. Then release the supports so that they stick to the foam insole. Slide the first strap from the top into the holes on both sides of the housing.
4. tighten all hooks and straps in order from the front of the foot to the top of the leg.

ADJUSTING AIR CHAMBER PRESSURE

- 1 Turn the button to the "pump" position, then depress the pump to inflate the chamber.
2. turn the button to the "let air out" position to deflate.

SIZE CHART (SHOE SIZE):

S: 36 - 39
M: 40 - 44
L: 44 - 46,5
XL: 46.5+

Length: AT53006 - 46 cm, AT53030 - 33 cm

SCOPE OF DELIVERY

Orthosis, instruction manual

NOTES:

In the event of a device-related "serious incident" that directly or indirectly led to, could have led to, or is likely to lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

NOTE: It is forbidden to use the product other than in accordance with its intended use

ATTENTION: The manufacturer is not responsible for damage caused by negligent maintenance, inadequate servicing or resulting from failure to follow the recommendations contained in this instruction manual.

TARGET PATIENT GROUP

The healthcare professional shall indicate, on his own responsibility, the use of the device for adults and children, taking into account the available variants/sizes/unnecessary functions/sizes and indications, bearing in mind the information provided by the manufacturer.

CLEANING AND MAINTENANCE

Remove metal parts and casing before washing (clean the casing with a damp cloth and mild detergent). Wash by hand in cold water using mild soap. Dry flat only in the shade. Do not tumble dry, dry clean or use bleach. Do not expose to sources of heat.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

HOW TO DISPOSE OF THE PRODUCT

The medical device can be disposed of as normal household waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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