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
AT03101, AT03101D, AT03102, AT03102D, AT03103, AT03103D

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
Anti-decubitus mattress

AT03101, AT03101D, AT03102, AT03102D, AT03103, AT03103D

Anti-decubitus foam mattresses for hospital beds. The mattresses allow optimal distribution of the pressure exerted by the body on the tissues adhering to the ground. Thanks to this, they slow down the process of clogging of blood vessels, leading to the formation of bedsores. They are recommended for people with severely limited mobility who spend most of their time lying down.

The appropriate elasticity of the foam from which the mattress is made ensures even distribution of body weight when lying down or changing position, preservation of the natural curvature of the spine, healthy sleep, prevention of spinal disorders.

Dimensions:

AT03101, AT03101D Foam anti-decubitus mattress (D-cover for disinfection): 200x90x10 cm AT03102, AT03102D Foam anti-decubitus mattress (D-cover for disinfection): 200x90x15 cm AT03103, AT03103D 7-STREF foam mattress (D-cover for disinfection): 200x90x15 cm

The maximum user weight is 130 kg.

Washing instructions for the waterproof (standard) cover:

1. Do not wash the foam!
2. Hand wash only the cover.
3. Clean using a damp sponge with laundry detergent.
4. Maximum washing temperature of 40°C.
5. Do not bleach! (Do not use free chlorine emitting compounds for bleaching).
6. Do not chemically clean!
7. Do not tumble dry.
8. Do not iron! (Do not allow contact with heated surfaces such as a radiator).
9. Deterioration of visual-aesthetic properties (surface wrinkling type) does not change the performance of the material (physical-mechanical properties).

Washing instructions for disinfecting cover:

1. Do not wash the foam!
2. Clean using a damp sponge with laundry detergent.
3. Maximum washing temperature of 95°C.
4. Do not bleach! (Do not use free chlorine emitting compounds for bleaching).
5. Do not chemically clean!
6. Do not tumble dry.
7. Do not iron! (Do not allow contact with heated surfaces such as a radiator).
8. Deterioration of visual-aesthetic properties (surface wrinkling type) does not change the performance of the material (physical-mechanical properties).

COMMENTS AND WARNINGS:

1. It is forbidden to use the product on children under 3 years old. The soft texture of the sponge can seriously interfere with the breathing process when the child lies in the stomach position.
2. Do not lie on the mattress in the stomach position.
3. Avoid contact of the product with fire. The product may catch fire.

ATTENTION: the mattress does not fully protect against the formation of bedsores in people who remain in a motionless lying position for a long period of time (e.g. paralyzed people, post-stroke, unconscious people, etc.). **NOTE:** it is forbidden to use the product other than in accordance with its intended use.

CONTRAINDICATIONS

physical or mental limitations (e.g., visual impairment) that prevent safe handling.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

NOTES:

The manufacturer is not responsible for damage caused by negligent cleaning or resulting from failure to follow the recommendations in this manual.

NOTES:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person or
- c) a serious threat to public health

the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions, or other disturbing symptoms that are unclear to you related to the use of the medical device, consult a healthcare professional.

NOTES:

The unfolded mattress takes about 72 h to fully achieve its properties, shape and full volume.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this device is intended (see the section on the purpose of the device of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Both in the case of purchasing the device independently, as well as the indication by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer, should be taken into account.

DISPOSAL OF THE PRODUCT

Dispose of the product in accordance with the applicable disposal conditions.

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER.

WARRANTY CARD

Model:

Serial Number:

Seal of the point of sale

and legible signature of Seller:

- 1) ANTAR Sp. J., 03-068 Warsaw, Zawisłańska 43 St., hereby provides a 12-month warranty for the product, starting from the date of purchase of the product by the Buyer. For products, the purchase of which is reimbursed by the National Health Fund, the warranty period is extended and covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
- 2) During the warranty period, ANTAR undertakes to perform the repairs necessary to restore the proper functioning of the product, within a period of 14 days from the date of receipt of the product by the service along with the proof of purchase and the warranty card.
- 3) ANTAR reserves the right to replace the product if it is determined that the repair costs are not cost-effective. Repair or replacement of defective products does not extend the warranty period.
- 4) If the product is found to be malfunctioning, it should be sent immediately to ANTAR's address (given in Section 1) or contact the point of sale.
- 5) Before use, read the instructions carefully and follow them. Use of the product contrary to its intended use and recommendations will result in irrevocable loss of warranty.
- 6) The warranty does not cover damage caused by external mechanical force, loading in excess of the user's maximum permissible weight as described in the manual, and contact of the product with alcohol, oily substances or gasoline. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
- 7) The product submitted to the service for a complaint must not cause an epidemic.
- 8) THIS WARRANTY IS VALID ONLY WITH PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A CLAIM, THE PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO RECOGNIZE THE EXTENDED WARRANTY PERIOD FOR PRODUCTS THE PURCHASE OF WHICH WAS REIMBURSED BY THE NATIONAL HEALTH FUND, A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES THAT ARE ORTHOPEDIC OBJECTS AND AIDS MUST BE ATTACHED.
- 9) Claims submitted without proof of purchase and warranty card with product serial number will not be accepted.


ANTAR Sp. J.
43 Zawisłańska Street
03-068 Warsaw, Poland
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net e-mail: antar@antar.net

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