AT51137
INSTRUKCJA UŻYTKOWNIA, CHODZIK
INSTRUCTION MANUAL, ROLLATOR
GEBRAUCHSANWEISUNG, ROLLATOR
NÁVOD K POUŽITÍ, CHODÍTKO
NÁVOD NA POUŽITIE, CHODÍTKO
MANUEL D'INSTRUCTION, DÉAMBULATEUR
GEbruIKSAANWIJZING, ROLLATOR
MANUAL DE INSTRUCCIONES, ANDADOR
MANUALE DI ISTRUZIONI, DEAMBULATORE
BRUKSANVISNING, RULLATOR

MD

PL

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EN

Thank you for purchasing our product. Please read the user manual carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, notes, and recommendations, contact a healthcare professional or the seller to avoid product damage or personal injury.

WARNING: Check all parts for any damage that may have occurred during transport. If such damage is noticed, DO NOT use the product. For more information, contact the seller.

TARGET PATIENT GROUP

Individuals suffering from illnesses, dysfunctions, or injuries for the treatment/rehabilitation or compensation of which this product is intended (see the section on product purpose in this manual). The product may be purchased independently by the user or based on recommendations from a doctor, therapist, or other specialist. Whether purchased independently or recommended by a healthcare professional, factors such as available sizes, necessary functions, product variants, indications and contraindications, as well as information provided by the manufacturer, should be taken into account.

APPLICATION / INDICATIONS

The product is designed to support mobility for individuals with limited walking ability. The rehabilitation rollator is designed for single-person use. It is suitable for:

Balance training: Allows patients to practice maintaining balance, which is crucial in rehabilitation after injuries or neurological disorders.

Posture support: Patients can train correct posture, especially important for individuals with postural defects or stability issues.

Improving motor coordination: Exercises performed with the rollator can help improve coordination and movement control, which is vital for patients with neurological conditions.

Rehabilitation after injuries: Assists in recovery following lower limb injuries, helping patients gradually regain mobility.

Increasing independence: For some patients, the rollator can serve as a tool to increase independence in daily life by enabling self-mobility.

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe use of the product.

Using a rehabilitation walker may be unsuitable for individuals with serious spinal stability problems or spinal injuries, as it may increase the risk of further injury.

In patients with uncontrolled seizures or epilepsy, use of the rollator may pose risks to both the patient and those nearby.

For individuals with very limited muscle strength, use of the rollator may be difficult or ineffective. In such cases, alternative rehabilitation methods, such as manual therapy or strength-supporting exercises, may be necessary.

For patients with severe balance disorders, using the rollator may increase the risk of falls and additional injuries.

COMPONENTS

Rollator, Bag, Backrest, User manual

PRODUCT FEATURES

Made of seamless tubing – lightweight, easy to store and transport.

Provides excellent maneuverability and stability.

Adjustable-height push handle – adaptable to the user.

Equipped with a shopping bag – allows carrying small items.

TECHNICAL SPECIFICATIONS

Frame: carbon fiber, water-transfer finish

Wheel fork: PA+TPE

Shopping bag: Oxford weave, 100% fabric, weight 0.28 kg

Seat: 100% Oxford fabric

Backrest: 100% Oxford fabric (same as seat)

Handles: PA+TPE, smooth non-slip ergonomic surface

Maximum load capacity: 150 kg

Maximum bag load: 5 kg

Total weight (without bag): 5.5 kg

Length: 65 cm

Total width: 65 cm

Total depth: 30 cm

Folded height (minimum): 84 cm

Unfolded width: 63 cm

Folded width: 22 cm

Folded dimensions (L × W × H): 65 × 23 × 85.5 cm

Seat width: 48 cm

Seat depth: 21 cm

Seat height: 62 cm

Distance between handles: 57 cm

Adjustable handle height (min–max): 76–100 cm

Brakes: 2 × maneuvering brakes with parking function

Front wheels: 8 in / 200 mm, solid TPR&PP tires, smooth tread, hub with sealed bearings

Rear wheels: 200 × 36 mm, black PA+15%GF+TPR, soft-surface tire, parking brake system (steel cable)


WARNING! Product dimensions may vary by ± 2 cm.

This symbol indicates the maximum user weight.

BACKREST ASSEMBLY

Attach the backrest by fastening the straps to the side tubes of the rollator.

UNFOLDING THE ROLLATOR

Release the flexible folding lock strap.

Hold the handles with both hands and pull them apart.

Press the support tubes under the seat cushion on both sides until you hear the locking mechanism click.

WARNING: After unfolding or assembling the rollator, always make sure that the seat support is securely locked and the frame does not move. Always ensure that both handle tubes are safely locked into the main support tube with the height adjustment locking button to prevent accidental sliding under pressure. Always check that all accessories are correctly and securely attached BEFORE use. Failure to follow these warnings may result in personal injury.

FOLDING THE ROLLATOR

Set the support rods on both sides to the lowest position.

Pull the folding lock strap located in the middle of the seat cushion until the left and right frames are fully folded.

Secure the seat support tubes on both sides with the elastic straps.

ADJUSTING HANDLE HEIGHT

Release the handle locking knob.


Set the height according to your needs. Both handles should be at the same height.

Properly adjusted handles should be at the user's hip level.

Never leave the device without engaging the brakes.

USING THE BRAKES

To stop the rollator, pull the brake levers upward toward the handles.

To remain stationary, push the brake levers downward into the “pressed” position.

To start moving, release the brake levers – they will automatically return, releasing the brakes.

Always apply the brakes when you do not intend to move the device.

Never leave the device without brakes applied when someone is sitting on it!

Never leave the device on a slope without brakes applied!


When sitting on the rollator, always sit upright, with forearms or hands lightly resting on the handles.

Never lean backward when sitting on the rollator.

Locking/Unlocking the brakes:

To lock the brakes, press the brake levers downward until you hear a “click.”

To unlock, pull the brake levers upward until you feel them loosen.

WARNINGS

If a serious incident related to the product occurs, which directly or indirectly results in, may result in, or could have resulted in any of the following:

- death of the patient, user, or another person,
- temporary or permanent deterioration of health of the patient, user, or another person,
- serious risk to public health,

this “serious incident” must be reported to the manufacturer and the competent authority of the member state where the user or patient resides. NOTE: In case of pain, allergic reactions, or any other concerning or unclear symptoms related to the use of the medical device, consult a healthcare professional.

NOTE: It is prohibited to use the product in any way other than its intended purpose.

LIMITATIONS OF USE

The seat must not be loaded with more than 150 kg.

The push handles must not be loaded with more than 150 kg.

The bag has a load limit of 5 kg.

The device is intended for use on flat surfaces, indoors and/or outdoors.

The device is not intended for self-propulsion in a seated position.

Brakes MUST be applied BEFORE sitting on the seat.

ACCIDENT PREVENTION

Use the device only as intended.

Use the device only in good technical condition.

Do not modify the construction of the device unless you have written permission from the manufacturer.

When making adjustments, be careful not to pinch fingers or body parts between moving parts.

Always use the rollator on a stable surface.

All wheels MUST remain in contact with the ground at all times to ensure stability and avoid accidents.

When using the rollator in a stationary position, the hand brakes MUST be applied.

CLEANING

Clean the device regularly and keep it tidy.

Light dirt and dust can be removed with a damp cloth.

Heavier dirt can be cleaned with a mild, non-corrosive detergent and warm water.

Do not use strong abrasives or chlorine-based cleaners as they may cause damage or discoloration.

Dry parts with a soft cloth.

Wheels: Clean wheels with warm water and mild detergent. Dry with a clean cloth. Lubricate periodically with multi-purpose oil or when wheels start to squeak.

MAINTENANCE

To ensure safe and proper use, inspect the device regularly according to the schedule below:

Weekly checks:

Check that all screws are tight.

Inspect brake cables for looseness or damage.

Tighten screws and nuts if necessary.

Monthly checks:

Check handle stability and condition.

Inspect brake levers for proper function.

Inspect bases (legs) for stability and looseness.

Quarterly checks (every 4 months):

Inspect wheels to ensure they rotate freely and are undamaged.

Check the bag for wear or damage.

Inspect the frame and supports for cracks or deformations.

TRANSPORT

The rollator folds easily, making it convenient to transport. For transportation, fold it as described above. Always transport the rollator securely to prevent sudden, unsafe movements.

REUSE

The rollator is suitable for reuse. Before passing it to another user, it must be inspected by a qualified dealer. Check the following:

Visual condition of the frame.

Secure fit of screws and nuts.

Condition of plastic parts.

Condition of wheels and bearings.

Function and condition of brakes.

If any part is damaged or broken, it must be replaced. Wheels should be replaced whenever significant wear is noticed. All loosened screws and nuts must be tightened. The product must be thoroughly cleaned and disinfected before reuse.

WARNING: Do not reuse the rollator if the frame is bent or broken.

The user manual must always accompany the rollator when changing users.

STORAGE

Store the product in a dry place. Do not expose it to direct sunlight.

DISPOSAL

When taken out of service, dispose of the product as regular household waste. Follow local waste separation regulations. Protect the environment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. AN TAR Sp. J., 03-068 Warsaw, ul. Zawiślńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to AN TAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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